

Leaf Spring Orthosis



Features and Benefits

- Heat mouldable splint providing dorsiflexion assistance and limited plantar flexion
- Easily trimmed to fit into a shoe and with padded touch and close fastening
- 3/4 and full foot versions are available

Indications for Use

- Flaccid foot drop

3/4 Foot

Size	Heel-Calf Height (cm)	Foot Length (cm)	UK Shoe Size	Product Code Left	Product Code Right
Small	33	15	2.5 - 5.5	PLSC/100/L	PLSC/100/R
Medium	35.5	15.5	5.5 - 8.5	PLSC/200/L	PLSC/200/R
Large	35.5	18.5	8.5 - 10.5	PLSC/300/L	PLSC/300/R
X-Large	38	29.5	10.5 - 12.5	PLSC/400/L	PLSC/400/R

Full Foot

Size	Heel-Calf Height (cm)	Foot Length (cm)	UK Shoe Size	Product Code Left	Product Code Right
Small	33	24	2.5 - 5.5	PLS/FF/100/L	PLS/FF/100/R
Medium	34.5	25.5	5.5 - 8.5	PLS/FF/200/L	PLS/FF/200/R
Large	37.5	28	8.5 - 10.5	PLS/FF/300/L	PLS/FF/300/R
X-Large	37.5	30.5	10.5 - 12.5	PLS/FF/400/L	PLS/FF/400/R

Dynamic AFO



Features and Benefits

The Dynamic Walk is a dynamic foot lifting AFO made of a lightweight carbon fibre composite and double sided peek rods that provide high strength and flexibility.

- This unique design allows patients to drive and climb stairs since it provides stable, free movement and constant lifting of the foot
- The open heel construction and slim foot design allows the AFO to fit easily into normal shoes
- Standard sizes below. Wide ankle also available, contact Customer Services for details

Indications for Use

- Flaccid foot drop

Size	Foot Length (cm)	Loaded Height (cm)	Calf Circ. (cm)	UK Shoe Size	Product Code Left Medial	Product Code Right Medial	Product Code Left Lateral	Product Code Right Lateral
Small	23	32	28.5	4	DW621	DW622	DW821	DW822
Medium	25	34	34.5	6	DW641	DW642	DW841	DW842
Large	27	37	39.5	8	DW661	DW662	DW861	DW862

Size	Foot Length (cm)	Loaded Height (cm)	Ankle Width (cm)	Calf Circ. (cm)	UK Shoe Size	Product Code Double Left	Product Code Double Right
Small	23	32	67.5	28.5	4	DW221	DW222
Medium	25	34	72.5	34.5	6	DW241	DW242
Large	27	37	80.5	39.5	8	DW261	DW262