

Sacro Lumbar Support



Features and Benefits

- Made from a ventilated mesh with four posterior shaped stays for greater support
- Elastic front with one-piece touch and close fastening; elasticated double fulcrum straps for additional support

Indications for Use

- Osteoarthritis
- Lower back pain
- Sciatica
- Post-op stability

Size	Waist Circ. (cm)	Product Code
Small	63 - 73	SLS/S
Medium	74 - 83	SLS/M
Large	84 - 96	SLS/L
X-Large	97 - 109	SLS/XL

Elasticated Lumbar Sacral Support



Features and Benefits

- Lumbar sacral support made from lightweight ventilated mesh elastic
- Two posterior shaped stays and two plastic stays
- Two touch and close fastenings allows the patient to obtain the optimum fit
- Elasticated double fulcrum straps for additional support and improved fit

Indications for Use

- Osteoporosis
- Lower back pain
- Sciatica

Size	Waist Circ. (cm)	Product Code
Small	69 - 76	ELS/S
Medium	79 - 86	ELS/M
Large	89 - 97	ELS/L
X-Large	99 - 107	ELS/XL

Depth - 24.1 cm

Deep Elasticated Lumbar Sacro Support



Features and Benefits

- Deep lumbar sacral support made from lightweight ventilated elastic
- Four posterior stays
- Adjustable touch and close fastening allows the patient to obtain optimum fit
- Elasticated double straps for additional support and improved fit

Indications for Use

- Osteoporosis
- Lower back pain
- Sciatica
- Osteoarthritis

Size	Waist Circ. (cm)	Product Code
Small	69 - 76	DEL/S
Medium	79 - 86	DEL/M
Large	89 - 97	DEL/L