

Universal PROGLIDE 174



Features and Benefits

The PROGLIDE 174 cervical collar is designed to offer maximum support with the addition of an extended sternal pad for increased occipital support. The easy glide system has 12 adjustments available, allowing for improved fit and aims to be a one-size-fits-all solution. The collar maintains posture and with improved padding, allows more effective pressure distribution to decrease the risk of redness and skin breakdown. The collar also combines a new and improved hook and loop fastening system, providing the clinician with advanced control whilst fitting.

- Extended sternal pad and occipital support
- Ergonomic adjustable control system
- 12 individual adjustments
- Pulley system cervical collar
- Removable thoracic attachment
- Soft cloth inner liner
- Tracheal opening
- MRI, CT, X-Ray compatible

Indications for Use

- Cervical spine precaution for trauma patients
- Immobilisation for pre and post cervical spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis
- Anterior cervical fusion
- Halo removal
- Dens type 1 cervical fractures of C2 vertebrae
- Anterior discectomy
- Teardrop fracture of the vertebral body
- Cervical strain
- Neural involvement
- Atlanto-occipital dislocation

Size	Description	Product Code
Universal	174 Collar	PG174
Universal	Spare liner set	PGPAD174

Miami J® Advanced



Features and Benefits

- Unique sizing mechanism allows for easy modification to suit individual patients, simultaneously achieving an accurate fit and reduced stock levels
- Bioengineered to minimise pressure points in key known areas: chin, occiput, trapezius and clavicle
- Sternal pad disperses ambulatory pressure and is fully removable, enabling the front of the collar to be taken off when cleaning intubated patients, without the need to remove the tube
- Patented Flex Edge® technology delivers soft over-moulded edges that gently conform to the contours of the anatomy and distribute pressure evenly
- MRI, CT and X-ray lucent


Indications for Use

- Post trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post operative

Phenotype Sizing Chart

Speciality sizes: 15% of the population

Standard sizes: 85% of the population



With Miami J® Advanced, one collar can be adjusted to fit various sizes, please contact Customer Services for guidance on sizing.

Description	Product Code - Collar Set*	Product Code - Collar
Miami J® Advanced	MJAR-101	MJA-100
Replacement pads	MJAP-101	

* The collar set comes with one additional set of replacement pads.