

OMO HIT



Features and Benefits

- Stabilisation of shoulder
- No stretching in shoulder-neck region
- Secure fixing of arm
- Early functional follow-up treatment eases post-operative physiotherapy
- Anatomically contoured support with silicone insert for friction massage on the musculofascial trigger points
- Cold/hot pad to reduce pain and oedema
- Axillary freedom allows easy body hygiene and increased comfort when wearing

Indications for Use

- Follow-up treatment after shoulder luxation, shoulder sub-luxation
- Sprain, contusion
- Fixing of various upper arm fractures
- Shoulder blade fracture
- Follow-up treatment in injuries to the rotator cuff
- Shoulder bursitis
- Post-operative treatment

Size	Chest Circ. (cm)	Product Code Left	Product Code Right
Small	< 92	07262/SL	07262/SR
Medium	92 - 102	07262/ML	07262/MR
Large	102 - 115	07262/LL	07262/LR
X-Large	115 +	07262/XLL	07262/XLR

Front Closure Clavicle Support



Features and Benefits

Padded straps cross in a figure of eight to provide secure fit and alignment of clavicle fractures.

Indications for Use

- Clavicle fracture
- Mild kyphosis

Size	Chest Circ. (cm)	Product Code
X-Small	51 - 61	205032
Small	61 - 76	205033
Medium	76 - 91	205035
Large	91 - 107	205037
X-Large	107 - 122	205038
XX-Large	122 - 137	205039

Subluxation Cuff



Features and Benefits

- Proven, effective design with comfortable padded straps
- Subluxation is resisted by the pull of the strapping acting on the upper arm across the shoulder
- Universal left/right application

Indications for Use

- Hemiplegia
- Shoulder subluxation

Size	Mid Upper Arm Circ. (cm)	Product Code
Small	22 - 27	SLC/S
Medium	27 - 31	SLC/M
Large	31 - 35	SLC/L
X-Large	35 - 40	SLC/XL