

VENUM Universal LSO



Features and Benefits

The VENUM LSO is the ultimate one-size-fits-all semi rigid comfort brace for the lumbar-sacral region. Its primary function is to reduce pain, whilst providing support and immobilisation to the lumbar-sacral spine. The highly effective cincher design produces maximal intra-abdominal pressure to distract the vertebral bodies of the spine, reducing pain, improving non/post-operative outcomes and enhancing patient rehabilitation. As the brace is universal in nature, it has the benefit of fitting every patient type and can be stored easily where limited storage is available.

- Double drawstring pulley closure system
- Cincher technology
- Reduced stock holding
- Cost effective
- Motion restriction of the lumbar regions L1 - L5
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic arthritic grips can be customised to patient size by extending or retracting cable length
- Made with durable, breathable, soft, hypo-allergenic materials that wick away perspiration
- Fits waist circumferences from 60cm-157.5cm; optional extender panel adds 10.25cm to circumference

Indications for Use

- Chronic back pain
- Degenerative, herniated or bulging disc
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- Spinal instability
- Burst Fx
- Pre or post-op spinal lumbar stabilisation
- Pre or post-op disk surgery

Size	Description	Product Code
Universal	VENUM LSO belt	VHBELT
Universal	Extender Panel	STEXT

Size Chart

To fit the VENUM LSO belt correctly, please choose the desired size on the slide panel from the charts below.

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5		

Oasis LSO



Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid LSO braces.

- The traction closure system provides easy compression adjustments
- Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates
- Flexibility of closing mechanism ensures an excellent fit across various body shapes
- Effective within T9 - S1

Indications for Use

- Burst Fx
- Post-op spinal lumbar
- Stabilisation
- Chronic back pain
- Compression Fx
- Degenerative, herniated or bulging disc
- Spinal instability
- Post-op disc surgery
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours or inflammation
- Post-op laminectomy

Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OALS/XS
Small	76.2 - 91.4	88.9 - 104.14	OALS/S
Medium	86.4 - 101.6	99 - 114.3	OALS/M
Large	96.5 - 111.8	109.2 - 124.5	OALS/L
X-Large	106.7 - 121.9	119.4 - 134.6	OALS/XL
XX-Large	116.8 - 132	129.5 - 144.8	OALS/XXL
XXX-Large	127 - 142.2	139.7 - 154.9	OALS/XXXL