

VENUM Universal TLSO PPX and VENUM Universal HYBRID TLSO Size Chart

To fit the VENUM TLSO PPX and the VENUM HYBRID TLSO correctly, please choose the desired size on the slide panel from the charts below. Note: the anterior side needs to be set at one size larger to allow for anterior attachment.

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5		

Oasis TLSO



Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid LSO braces. Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates.

The Oasis TLSO and TLSO LP 3 Panel DLK includes a Dorsal Lumbar Kit (DLK) with Sternal Panel to prevent forward flexion. Effective within T4 - S1.

Indications for Use

- Burst Fx
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Kyphosis
- Laminectomy
- Osteoporosis
- Paraplegia
- Post-operative
- Spinal stabilisation
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis

Model wearing SPK version shown adjacent. (SPX ATE version fitted to the Edge model shown on next page)



Size	Waist (cm)	Hip (cm)	Product Code	Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OATL4-XS-PPK	X-Small	66 - 82.3	78.7 - 94	OATL4-XS-SPK
Small	76.2 - 91.4	88.9 - 104.14	OATL4-S-PPK	Small	76.2 - 91.4	88.9 - 104.14	OATL4-S-SPK
Medium	86.4 - 101.6	99 - 114.3	OATL4-M-PPK	Medium	86.4 - 101.6	99 - 114.3	OATL4-M-SPK
Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-PPK	Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-SPK
X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-PPK	X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-SPK
XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-PPK	XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-SPK
XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-PPK	XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-SPK

Available in SPK & PPK versions, other versions available to order.