

## UFO™



### Features and Benefits

Introducing an improved prefabricated night splint for the effective treatment of plantar fasciitis and achilles tendinitis. The UFO™ comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Load Prolonged Stretch [LLPS] to the soft tissues of the plantar fascia and achilles tendon. LLPS promotes gentle stretching, healing and re-modeling of these tissues, effectively reducing the stress and pain caused by inflammation and tissue shortening. Such support can help break the cycle of injury, pain and re-injury due to otherwise unlimited contraction and deformity of the soft tissues at night.

- The simple variable dorsiflexion design allows the patient to quickly and easily adjust the amount of applied stretch to provide the most effective pain relief
- Adjustable bilateral dorsiflexion assist straps
- CoolFoam™ replaceable liner wicks moisture away from the skin for enhanced hygiene

### Indications for Use

- Plantar fasciitis
- Achilles tendonitis

| Size   | Footplate Length (cm) | Brace Height (cm) | Product Code |
|--------|-----------------------|-------------------|--------------|
| Small  | 14                    | 16.5              | B/3614       |
| Medium | 17                    | 25                | B/3615       |
| Large  | 19                    | 26.5              | B/3616       |
| Custom | N/A                   | N/A               | B/310/PE     |

## Paediatric A-Boot



### Features and Benefits

- Washable, soft fleece-lined cover wicks away moisture away from the skin to maintain skin integrity and allows heel to float, thus eliminating pressure of friction on the heel and enhancing blood circulation
- Rubber sole with non-slip surface provides a walking base so patients can walk comfortably
- The splint frame resists plantar flexion while providing a dynamic push toe extension for the correction of foot drop, foot ankle contractures and deformity
- The adjustable toe post relieves pressure on the toes and adjusts to accommodate a variation of foot sizes
- Rotation bar can be positioned to the side to control hip and leg rotation
- An extender kit is available for patients with increased calf and foot circumference and or edema
- Universal sizing to reduce stock holding

### Indications for Use

- For the patient who requires positioning but can also transfer and ambulate

| Size      | Foot Length (cm) | Calf Circ. (cm) | Heel to Top of Back (cm) | Product Code |
|-----------|------------------|-----------------|--------------------------|--------------|
| Universal | 12.7 - 20.3      | Up to 17.9      | 17.78                    | A-BOOT-P     |

## Universal PRAFO® Paediatric



### Features and Benefits

The PRAFO® (Pressure Relief Ankle Foot Orthosis) is a totally adjustable, custom fitted AFO that can help manage many ankle/foot anomalies, offering a range of adjustable features.

- Custom-contoured aluminium heel connector bar helps to control dorsi-plantar flexion by allowing infinite, measurable adjustments
- Ergonomically designed foot plate allows the liner to be changed easily without tools
- Integral, securely fastened vibram walking base provides a shoe-type surface for ambulatory patients
- The design protects the patient around the vulnerable areas of the foot and heel and will remain so during ambulation
- Calf and foot segments can be custom-modified to accommodate patient physiology
- Anti-rotation bar made of malleable aluminium and adjusts to control unwanted extremity rotation
- Suitable for day and night use
- Universal sizing to reduce stock holding

### Indications for Use

- Pressure ulcers
- Can be used to resist the development of contracture

| Size      | Product Code |
|-----------|--------------|
| Universal | PRFO-550     |