

## Philadelphia Collar



### Features and Benefits

- Two-piece moulded collar provides adjustable and comfortable support
- Touch and close fastening


**Material:** Non-toxic hypo-allergenic plastazote.

### Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis
- Post-op control of extension/flexion

		Product Code	
Size	Neck Circ. (cm)	Front Depth 8.5 cm	Front Depth 11 cm
Small	25 - 33	PCC/8.5/S	PCC/11/S
Medium	33 - 41	PCC/8.5/M	PCC/11/M
Large	41 - 49	PCC/8.5/L	PCC/11/L

## Universal PROGLIDE 172



### Features and Benefits

The PROGLIDE 172 cervical collar is one of our most advanced collars to date. It maintains posture and with improved padding, allows more effective pressure distribution to decrease the risk of redness and skin breakdown. The easy glide system has 12 adjustments available, allowing for improved fit and aims to be a one-size-fits-all solution. The collar also combines a new and improved hook and loop fastening system, providing the clinician with advanced control whilst fitting.

- Ergonomic adjustable control system
- 12 individual adjustments
- Pulley system cervical collar
- Removable thoracic attachment
- Soft cloth inner liner
- Tracheal opening
- MRI, CT, X-Ray compatible



### Indications for Use

- Cervical spine precaution for trauma patients
- Immobilisation for pre and post cervical spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis
- Anterior cervical fusion
- Halo removal
- Dens type 1 cervical fractures of C2 vertebrae
- Anterior discectomy
- Teardrop fracture of the vertebral body
- Cervical strain
- Neural involvement
- Atlanto-occipital dislocation

Size	Description	Product Code
Universal	172 Collar	PG172
Universal	Spare liner set	PGPAD172