

## Kraemer Heel Cups



### Features and Benefits

- Children's multi-layered prefabricated foot orthosis covered in Multiform A30 shore
- Thermo-mouldable reinforced base and wedge

### Indications for Use

- Ligament weakness
- Sprains/strains
- Planovalgus feet

Shoe Size (EU)	Product Code	Shoe Size (EU)	Product Code
19 - 20	ZKR1-0607-19/20	31 - 32	ZKR1-0640-31/32
21 - 22	ZKR1-0607-21/22	33 - 34	ZKR1-0640-33/34
23 - 24	ZKR1-0607-23/24	35 - 36	ZKR1-0640-35/36
25 - 26	ZKR1-0607-25/26	37 - 38	ZKR1-0640-37/38
27 - 28	ZKR1-0607-27/28	39 - 40	ZKR1-0640-39/40
29 - 30	ZKR1-0640-29/30	41 - 42	ZKR1-0640-41/42

## EZI-KID (Peddy) Insoles



### Features and Benefits

- Dynamic insole for children
- Modular sandwich-like construction and flexible thermoformable wedge
- Soft sole to support the natural movement of the foot

### Indications for Use

- Back pain
- Posture related symptoms
- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Planovalgus feet

Shoe Size (EU)	Product Code	Shoe Size (EU)	Product Code
19 - 20	ZKR1-0640-19/20	29 - 30	ZKR1-0640-29/30
21 - 22	ZKR1-0640-21/22	31 - 32	ZKR1-0640-31/32
23 - 24	ZKR1-0640-23/24	33 - 34	ZKR1-0640-33/34
25 - 26	ZKR1-0640-25/26	35 - 36	ZKR1-0640-35/36
27 - 28	ZKR1-0640-27/28	37 - 38	ZKR1-0640-37/38
		39 - 40	ZKR1-0640-39/40

## Prostep



### Features and Benefits

Prostep is designed to restore the natural position of the talus by supporting the inside edge of the heel and supports the spring ligament. This resists the valgus rotation of the rearfoot, which is associated with pronation.

- Anti-bacterial top cover and forefoot and rearfoot integral posts of 0°, 3° & 5°
- Full length or 3/4 length available

### Indications for Use

- Low arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	3/4 Medium	3/4 Firm	Full Medium	Full Firm
Kids Small	13.5 - 1.5	PSSLM-KS	PSSLF-KS	PSFLM-KS	PSFLF-KS
Kids Large	2 - 3	PSSLM-KL	PSSLF-KL	PSFLM-KL	PSFLF-KL
X-Small	3.5 - 4.5	PSSLM-XS	PSSLF-XS	PSFLM-XS	PSFLF-XS
Small	5 - 6	PSSLM-S	PSSLF-S	PSFLM-S	PSFLF-S
Medium	6.5 - 7.5	PSSLM-M	PSSLF-M	PSFLM-M	PSFLF-M
Large	8 - 9	PSSLM-L	PSSLF-L	PSFLM-L	PSFLF-L
X-Large	9.5 - 10.5	PSSLM-XL	PSSLF-XL	PSFLM-XL	PSFLF-XL
XX-Large	11 - 12	PSSLM-XXL	PSSLF-XXL	PSFLM-XXL	PSFLF-XXL

## Kids Feet In Motion



### Features and Benefits

These brightly coloured pairs consist of one blue (left) and one right (red) to help children learn whilst providing a functional solution.

- Medium density EVA orthosis specifically designed for children's feet
- A flat base on the standard orthotic allows for modification with wedges and components
- Insoles are available as standard or with a pre-posted heel to 5°

### Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Planovalgus feet
- Postural related symptoms

UK Shoe Size	EU Shoe Size	Standard Product Code	3° Product Code	5° Product Code
4	19 - 20	SC3304	SC3330	SC3350
5	21 - 22	SC3305	SC3331	SC3351
6.5	23 - 24	SC3306	SC3332	SC3352
8	25 - 26	SC3307	SC3333	SC3353
9 - 10	27 - 28	SC3308	SC3334	SC3354
10 - 12	29 - 30	SC3309	SC3335	SC3355
13 - 1	31 - 32	SC3310	SC3336	SC3356
1 - 2	33 - 34	SC3311	SC3337	SC3357

## Trio



### Features and Benefits

The Trio orthotic is designed to resist excessive valgus rotation, by inverting the calcaneus within the heel cup thereby maintaining mid tarsal joint. It also supports the proximal medial arch of the foot providing navicular stability.

- Each pair is supplied with additional 3° medial biplanar rearfoot posts and 3° lateral forefoot
- Available in 3/4 and full length posts

### Indications for Use

- Plantar fasciitis
- Heel pain
- Calf strain
- Antero-medial knee pain
- Hip instability

Size	UK Shoe Size	3/4 Length Medium Density	3/4 Length Firm Density	3/4 Length Extra Firm Density	Full Length Medium Density	Full Length Firm Density	Full Length Extra Firm Density
Kids Small	12.5 - 13.5	TSLM-KS	TSLF-KS	TSLE-KS	TFLM-KS	TFLF-KS	TFLE-KS
Kids Large	1 - 2.5	TSLM-KL	TSLF-KL	TSLE-KL	TFLM-KL	TFLF-KL	TFLE-KL
X-Small	3 - 4.5	TSLM-XS	TSLF-XS	TSLE-XS	TFLM-XS	TFLF-XS	TFLE-XS
Small	5 - 6.5	TSLM-S	TSLF-S	TSLE-S	TFLM-S	TFLF-S	TFLE-S
Medium	7 - 8.5	TSLM-M	TSLF-M	TSLE-M	TFLM-M	TFLF-M	TFLE-M
Large	9 - 10.5	TSLM-L	TSLF-L	TSLE-L	TFLM-L	TFLF-L	TFLE-L
X-Large	11 - 12.5	TSLM-XL	TSLF-XL	TSLE-XL	TFLM-XL	TFLF-XL	TFLE-XL
XX-Large	13 - 14.5	TSLM-XXL	TSLF-XXL	TSLE-XXL	TFLM-XXL	TFLF-XXL	TFLE-XXL

## Trio Easy Fit



### Features and Benefits

Trio Easy Fit is intrinsically a 3/4 length medially posted orthotic with an additional posting option for triplanar control of pronation featuring a lateral cutaway design. This model maintains the essential features for pronation control, including an optional biplanar rearfoot post. A medial forefoot post can be used to replace the first metatarsal cushion if required, to achieve comfort without compromising on natural movement.

### Indications for Use

- Medium arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density
Kids	1 - 2.5	TEF-K
X-Small	3 - 4.5	TEF-XS
Small	5 - 6.5	TEF-S
Medium	7 - 8.5	TEF-M
Large	9 - 10.5	TEF-L
X-Large	11 - 12.5	TEF-XL

## Starflex Blanks



## Features and Benefits

A versatile polyester resin orthotic blank that is worked with hot water or a hot air gun.

- Available pre-cut to shoe size and can be moulded to a customised shape
- Very strong, yet thin and lightweight foot orthotic
- A 'posting' can be applied simply using adhesive

## Indications for Use

- Back pain
- Posture related symptoms
- Arch pain
- Knee pain
- Tendon pain

Paediatric Insole Sizes		Adult Sizes			
UK Shoe Size	Product Code	UK Shoe Size	Product Code	UK Shoe Size	Product Code
6.5	SC2014	3.5	SC2028	9	SC2042
8	SC2016	4.5	SC2030	9.5	SC2044
10	SC2018	5	SC2032	10.5	SC2046
11.5	SC2020	5.5	SC2034	11	SC2048
13	SC2022	6	SC2036	11.5	SC2050
1.5	SC2024	7	SC2038	12	SC2052
2.5	SC2026	8	SC2040		

## Peapod Junior



## Features and Benefits

- 1.5mm thermoplastic base offers excellent fitting compatibility in most footwear
- Heat mouldable and grindable

## Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code	UK Shoe Size	Product Code
5	SC2502	11.5	SC2510
6.5	SC2504	13	SC2512
8	SC2506	1.5	SC2514
10	SC2508	3.5	SC2516

## Control Tech Soft 3/4 Kids



## Features and Benefits

- Rearfoot wedging reduces the need for excessive arch height whilst maintaining support, comfort and providing superior shock absorption

## Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code
9 - 10	SC3051
10 - 12	SC3052
13 - 1	SC3053
1 - 2	SC3054

## Neurotec



### Features and Benefits

- Neurotec is designed to provide support under the medial arch and the peroneal notch, which gives stability to the subtalar joint and de-stabilises inversion and eversion
- The metatarsal head arch stabilises the forefoot and offloads the metatarsal heads in the forefoot
- The heel cup stabilises the heel and the toe crest dorsiflexes the toes to improve tone

### Indications for Use

- Mid neurological deficit affecting the foot
- Habitual reflex deformities
- Hypertonic muscle imbalance

Size	EU Shoe Size (Shell Length cm)	Product Code	Size	EU Shoe Size (Shell Length cm)	Product Code
A	11 (7.5)	NEUROA	K	21 (13.8)	NEUROK
B	12 (8.1)	NEUROB	L	22 (14.4)	NEUROL
C	13 (8.7)	NEUROC	M	23 (15.1)	NEUROM
D	14 (9.3)	NEUROD	N	24 (15.7)	NEURON
E	15 (10)	NEUROE	O	25 (16.3)	NEUROO
F	16 (10.6)	NEUROF	P	26 (16.9)	NEUROP
G	17 (11.2)	NEUROG	Q	27 (17.7)	NEUROQ
H	18 (11.8)	NEUROH	R	28 (18.3)	NEUROR
I	19 (12.6)	NEUROI	S	29 (18.9)	NEUROS
J	20 (13.2)	NEUROJ	T	30 (19.5)	NEUROT

## Podo Tech Junior - 3/4 Length



### Features and Benefits

- Designed especially for children
- Easily ground and heat mouldable

### Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size (EU)	Product Code	UK Shoe Size (EU)	Product Code
4 (19 - 20)	SC3220	13 (31 - 32)	SC3232
5 (21 - 22)	SC3222	1.5 (33 - 34)	SC3234
6.5 (23 - 24)	SC3224	3.5 (35 - 36)	SC3236
8 (25 - 26)	SC3226	5 (37 - 38)	SC3238
10 (27 - 28)	SC3228	6 (39 - 40)	SC3240
11.5 (29 - 30)	SC3230	8 (41 - 42)	SC3242