

Open Toe Cast Sandal



Features and Benefits

- Features open toes and heel for comfort and breathability
- Constructed from ribbed rocker sole with durable upper and contact closure

Indications for Use

- Lower extremity casted fractures

Size	Foot Length (cm)	Product Code
XXX-Small	10.2	69-3078
XX-Small	12.7	69-3079
X-Small	17.8	69-3080

Universal Rigid Post-Op Shoe



Features and Benefits

- Features flat footbed with rigid sole to prevent excessive foot motion whilst walking
- Constructed of wood-like sole with water resistant upper and contact closure
- Universal sizing to reduce stock holding

Indications for Use

- Post-operative foot protection and support

Size	Foot Length (cm)	Product Code
Universal	20.32	69-3121

Paediatric Fixed Walker



Features and Benefits

- Easy to apply design provides stability for the lower leg with low natural heel height
- Removable washable foam liner and aluminium uprights for adjustability; fits children from 1½ to 9 years old
- Constructed of coated malleable aluminium uprights, padded liner and contact closure straps

Indications for Use

- Stable foot or ankle fractures
- Severe sprains or strains
- Post cast removable or as a comfortable alternative to casting

Size	Footplate Length (cm)	Age (Years)	Product Code
Small	> 15.2	1 to 2	79-1201
Large	15.2 - 19	2 1/2 to 6	79-1202
X-Large	19 - 21.6	6 1/2 to 9 1/2	79-1203

