

A-Boot Universal



Features and Benefits

- Ambulation version (A-Boot) comes with removable, soft, non-slip sole
- Antimicrobial surface protector to shield against mould and mildew, stains, bacteria, fungus and algae
- Super soft fleece lined headliner cover is removable and machine washable
- A toe extension is available to accommodate all foot sizes
- Rotator bar controls hip and leg rotation
- Boots position heel in floating position to prevent pressure areas through pressure reduction
- Comes with extender kit to accommodate for wider or oedematous feet
- Universal sizing to reduce stock holding

Indications for Use

- Wider or oedematous feet
- Ankle contractures
- Heel ulcerations

Size	Product Code
Universal	A-BOOT

Also available in paediatric sizing - page 92

Open Toe Cast Sandal



Features and Benefits

- Open toes and heel for comfort and breathability
- Constructed from ribbed rocker sole with durable upper and contact closure

Indications for Use

- Lower extremity casted fractures

Size	Foot Length (cm)	Product Code
Small	23	69-3081
Medium	28	69-3082
Large	31	69-3083
X-Large	36	69-3084

All Purpose Boot



Features and Benefits

- Removable EVA insole is fully customisable to offload pressure points
- Wider sole and deep opening provides increased room to accommodate large dressings, bandages and casts
- Wide ankle strap eliminates heel slippage as well as securing the foot or cast in the shoe
- A rocker bottom sole encourages smooth roll through gait pattern and reduces plantar pressure

Indications for Use

- Post-operative protection for soft tissue trauma or surgery
- Plantar ulcers/wounds
- Closed toe allows dressing to keep clean and dry
- Wider sole allows more room for dressings and casts

Size	UK Shoe Size	Product Code
Small	4 - 5.5	APB-S
Medium	6 - 7.5	APB-M
Large	7.5 - 8.5	APB-L
X-Large	9 +	APB-XL

Post-Op Shoe

1



Features and Benefits

- Durable upper with soft padded lining
- Grooved non-slip flexible sole
- Touch and close through loop fastening
- Universal left/right application

Indications for Use

- Post surgery accommodation of bulky bandages

2

3

Size	UK Shoe Size	Product Code
Small	6 - 7	POS/S
Medium	7 - 9	POS/M
Large	9 - 10	POS/L
X-Large	10 +	POS/XL

4

Orthowedge Shoe



Features and Benefits

- Promotes healing by reducing weight from the forefoot
- Hard wearing padded nylon upper
- Adjustable touch and close fastening
- Universal left/right application

Indications for Use

- Ulcerations of the forefoot
- Painful forefoot
- Post surgery

5

6

Size	UK Shoe Size	Product Code
Small	3 - 5	OWS/S
Medium	5 - 7	OWS/M
Large	7 - 9	OWS/L
X-Large	9+	OWS/XL

7

Heel Wedge Heeling Shoe



Features and Benefits

- Forefoot closure provides uniform security
- Hard wearing padded nylon upper
- Adjustable touch and close fastening
- Universal left/right application

Indications for Use

- Ulcerations of the hindfoot

8

9

Size	UK Shoe Size	Product Code
Small	3 - 5	HWHS/S
Medium	5 - 7	HWHS/M
Large	7 - 9	HWHS/L
X-Large	9 +	HWHS/XL

Pressure Shoe



Features and Benefits

- Can be used as a post operative cast shoe, providing more protection for post operative procedures and fiberglass casts
- Touch and close fastening
- Extra double touch and close strap for security
- Durable sole material
- Universal left/right application

Indications for Use

- Post-op cast shoe
- Accommodation of dressings

Size	UK Shoe Size	Product Code
X-Small	> 4	PS/XS
Small	4 - 5.5	PS/S
Medium	6 - 7.5	PS/M
Large	7.5 - 8.5	PS/L
X-Large	9 +	PS/XL

Neoprene Boot



Features and Benefits

- 4mm neoprene rubber faced and lined with jersey nylon
- Touch and close fastening
- EVA unit with 8mm sole and 25mm heel
- Comes with heel stiffeners
- Comfort foam insole
- Available in black in wide fitting in full sizes between 2 - 12, in pairs only

Indications for Use

- Bulky dressing and fluctuating oedema

UK Shoe Size	Product Code	UK Shoe Size	Product Code	UK Shoe Size	Product Code
2	SHNB/02	6	SHNB/06	10	SHNB/10
3	SHNB/03	7	SHNB/07	11	SHNB/11
4	SHNB/04	8	SHNB/08	12	SHNB/12
5	SHNB/05	9	SHNB/09	13	SHNB/13

Felt Boot



Features and Benefits

- Supplied with heel and toe stiffeners
- Touch and close fastening
- EVA unit with 8mm sole and 25mm heel
- Comfort foam insole
- Available in black in wide fitting in full sizes between 2 - 12 in pairs only

Indications for Use

- Bulky dressing and fluctuating oedema

UK Shoe Size	Product Code	UK Shoe Size	Product Code	UK Shoe Size	Product Code
2	SHFB/02	6	SHFB/06	10	SHFB/10
3	SHFB/03	7	SHFB/07	11	SHFB/11
4	SHFB/04	8	SHFB/08	12	SHFB/12
5	SHFB/05	9	SHFB/09	13	SHFB/13