

## Hinged Elbow Sleeve



### Features and Benefits

- Metal hinges limit medial and lateral movement and restrict hyperextension
- Universal left/right application

Material: Coolmax®

### Indications for Use

- Bursitis
- Osteoarthritis
- Tendonitis
- Rheumatoid arthritis
- Joint inflammation

Size	Elbow Joint Circ. (cm)	Product Code
Small	23 - 25.5	HES/BX/BE/S
Medium	25.5 - 28	HES/BX/BE/M
Large	28 - 30.5	HES/BX/BE/L
X-Large	30.5 - 33	HES/BX/BE/XL

## ROM Elbow Brace



### Features and Benefits

- Soft foam wraps
- ROM hinge for setting flexion/extension
- Touch and close fastening, also supplied with soft padded arm sling

### Indications for Use

- Elbow contracture management
- Post-op elbow flexion/extension control

Size	Elbow Joint Circ. (cm)	Length (cm)	Product Code
Small/Medium	23 - 31	39	REB/S-M
Large/X-Large	32 - 40	39	REB/L-XL

1

2

3

4

5

6

7

8

9