

Standard Universal ROM Knee Brace



Features and Benefits

- Open four panel design
- Easily adjustable ROM hinges for setting flexion/extension
- Touch and close fastening
- Universal sizing to reduce stock holding

Indications for Use

- Post-op immobilisation of the knee
- Unstable knee
- Limitation of knee movement

Size	Product Code
Universal	ROM/L/U

Universal ROM Knee Brace



Features and Benefits

- The Universal ROM Knee Brace features butterfly padding to prevent migration
- With an adjustable length from 57-69cm, the fully adjustable Range of Motion (ROM) knee mechanism has a Drop-Lock system that can fix the joint at 0°, 15°, 20°, 30° and 45°
- Universal sizing to reduce stock holding
- Easy 'Slide & Set' ROM settings:
 - Flexion settings: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
 - Extension settings: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Compliance ROM protection – easy-fit plastic ties supplied to prevent altering prescribed ROM settings
- Telescopic push-button adjustment mechanism with ability to reposition above and below knee pads
- Full adjustable soft circumferential straps with return buckles and non-slip foam wraps for excellent patient comfort and fit

Indications for Use

- Improved joint function
- Post-operative knee
- Limitation of knee movement

Size	Product Code
Universal	57260

Contender Universal Post-op Knee Brace



Features and Benefits

- The Contender™ Post-op Knee Brace features a simple adjustment system to minimise fit time after surgery, and at 765g, it is one of the lightest braces on the market.
- The X-Flex™ Cuff System flexes to hug the leg, significantly reducing brace migration. QuikZip™ buckles allow for easy application and micro adjustment of the strap for a secure fit
 - The ZipTrak Telescoping System adjusts brace length from 45.7cm - 66cm, in 3mm increments and also allows adjustment of intermediate cuffs to help avoid incision sites
 - The hinge features a six position Drop-Lock that fixes the brace from -10° to 40° in 10° increments. Controlled range of motion is available from -10° to 110° and is easily set using flexion and extension finger touch buttons
 - Offers the latest in knee brace technology such as multi-telescoping cuffs, user friendly buckles and critical suspension systems
 - Length adjustable 45.7cm - 66cm. Fits up to 86.4cm thigh circumference. Measurement taken from leg circumference 15.2cm above mid patella

Indications for Use

- ACL, PCL, MCL and LCL surgeries
- Meniscal repairs
- Patella realignment
- Regenerative chondroplasty
- Stable femoral fractures
- Total knee replacements
- Tibial plateau fractures

Size	Product Code
Universal	75-7500