

Steeper motionSupport® Insoles



Features and Benefits

- motionSupport® insoles are a durable, innovative range of ready-made insoles featuring unique design and construction
- The full length polyurethane insoles offer improved forefoot control whilst the 3/4 length fibreglass reinforced shell and fibreglass extension on the lateral aspect permits a more propulsive gait
- The flexible construction under the first metatarsal head allows the first ray to plantarflex and a small metatarsal button reduces pressure on the metatarsal heads
- Available in three arch heights to optimise sagittal plane control
- The slimmer profile allows for easier fitting in footwear

Indications for Use

- Arch pain
- Over pronation
- Posture related symptoms
- Back pain
- Tendon pain
- Knee pain
- Shin splints
- Metatarsalgia
- Function hallux limitus

EU Shoe Size	Low Arch (Brown)	Normal Arch (Grey)	High Arch (Blue)
36 - 37	KRI-7055-36-37	KRI-7052-36-37	KRI-7053-36-37
38 - 39	KRI-7055-38-39	KRI-7052-38-39	KRI-7053-38-39
40 - 41	KRI-7055-40-41	KRI-7052-40-41	KRI-7053-40-41
42 - 43	KRI-7055-42-43	KRI-7052-42-43	KRI-7053-42-43
44 - 45	KRI-7055-44-45	KRI-7052-44-45	KRI-7053-44-45
46 - 47	KRI-7055-46-47	KRI-7052-46-47	KRI-7053-46-47

Steeper Support Hallux Rigidus Insoles



Features and Benefits

- Hallux Rigidus insoles are designed to treat painful extension at the first metatarsophalangeal joint
- An ortholyn inner shell extends distally beyond the first metatarsophalangeal joint providing a responsive shock absorbing core whilst limiting painful extension motion
- Choice of three arch profiles optimises the contouring to the foot
- The low profile and durable construction facilitates use in a wide range of footwear
- Insole design allows use with various heel heights
- Sizes outside of the chart below are available to order on request

Indications for Use

- Painful extension at the first metatarsophalangeal joint

Size	Product Code - Male			Product Code - Female		
	Low Arch	Normal Arch	High Arch	Low Arch	Normal Arch	High Arch
36	-	-	-	KRI-7405-36	KRI-7402-36	KRI-7403-36
37	-	-	-	KRI-7405-37	KRI-7402-37	KRI-7403-37
38	-	-	-	KRI-7405-38	KRI-7402-38	KRI-7403-38
39	-	-	-	KRI-7405-39	KRI-7402-39	KRI-7403-39
40	KRI-7415-40	KRI-7412-40	KRI-7413-40	KRI-7405-40	KRI-7402-40	KRI-7403-40
41	KRI-7415-41	KRI-7412-41	KRI-7413-41	KRI-7405-41	KRI-7402-41	KRI-7403-41
42	KRI-7415-42	KRI-7412-42	KRI-7413-42	KRI-7405-42	KRI-7402-42	KRI-7403-42
43	KRI-7415-43	KRI-7412-43	KRI-7413-43	-	-	-
44	KRI-7415-44	KRI-7412-44	KRI-7413-44	-	-	-
45	KRI-7415-45	KRI-7412-45	KRI-7413-45	-	-	-
46	KRI-7415-46	KRI-7412-46	KRI-7413-46	-	-	-
47	KRI-7415-47	KRI-7412-47	KRI-7013-47	-	-	-

Steeper Support Insoles



Features and Benefits

- Steeper Support Insoles feature an extremely low profile, durable construction with three arch profiles optimising the contouring to the foot
- An ortholyn inner shell provides a responsive shock-absorbing core whilst the metatarsal button relieves pressure on the metatarsal heads
- The option to select the male or female profile further optimises the individual fit
- Insole design allows use with variable heel heights
- Sizes outside of the chart below are available to order on request

Indications for Use

- Morton's neuroma
- Metatarsalgia
- Forefoot pain
- Mechanical foot pain

Size	Product Code Male			Product Code Female		
	Low Arch	Normal Arch	High Arch	Low Arch	Normal Arch	High Arch
36	-	-	-	KRI-7005-36	KRI-7002-36	KRI-7003-36
37	-	-	-	KRI-7005-37	KRI-7002-37	KRI-7003-37
38	-	-	-	KRI-7005-38	KRI-7002-38	KRI-7003-38
39	-	-	-	KRI-7005-39	KRI-7002-39	KRI-7003-39
40	KRI-7015-40	KRI-7012-40	KRI-7013-40	KRI-7005-40	KRI-7002-40	KRI-7003-40
41	KRI-7015-41	KRI-7012-41	KRI-7013-41	KRI-7005-41	KRI-7002-41	KRI-7003-41
42	KRI-7015-42	KRI-7012-42	KRI-7013-42	KRI-7005-42	KRI-7002-42	KRI-7003-42
43	KRI-7015-43	KRI-7012-43	KRI-7013-43	-	-	-
44	KRI-7015-44	KRI-7012-44	KRI-7013-44	-	-	-
45	KRI-7015-45	KRI-7012-45	KRI-7013-45	-	-	-
46	KRI-7015-46	KRI-7012-46	KRI-7013-46	-	-	-
47	KRI-7015-47	KRI-7012-47	KRI-7013-47	-	-	-

Diaped Duosoft Plus



Features and Benefits

- Therapeutic insoles for insensate feet
- The dual-layer design of Duosoft Plus consists of a resilient, shock absorbing urethane base and a PORON 94[®] top cover
- The slow rebound quality of the material cradles the foot and allows enhanced redistribution of body weight over the foot's plantar surface
- Soft, comfortable and allows the feet to remain cool, dry and odourless
- A Vygel heel accompanies the dual-layer formation to provide optimal shock attenuation at heel strike

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- Skin breakdown

UK Shoe Size	EU Shoe Size	Product Code
4 - 5	37 - 39	SB2402
6 - 7	40 - 41	SB2404
8 - 9	42 - 43	SB2406
9 - 10	44 - 45	SB2408
11 - 12	46 - 48	SB2410