

Long Wrist Brace



Features and Benefits


- Longer two-piece design giving support and comfort for continuous use
- Shape allows increased finger mobility and minimal thumb web discomfort
- Fully bound edges

Indications for Use

- Flaccid drop wrist
- Rheumatoid arthritis
- Ulnar and radial stability

Material: Coolmax®

Size	Wrist Circ. (cm)	Product Code - Coolmax® Beige	
		Left	Right
X-Small	Up to 14	LWB/CX/BE/XS/L	LWB/CX/BE/XS/R
Small	14 - 16	LWB/CX/BE/S/L	LWB/CX/BE/S/R
Medium	16.5 - 19	LWB/CX/BE/M/L	LWB/CX/BE/M/R
Large	19 - 21.5	LWB/CX/BE/L/L	LWB/CX/BE/L/R
X-Large	Over 21.5	LWB/CX/BE/XL/L	LWB/CX/BE/XL/R



About Coolmax®

CoolMax is the industry leader when it comes to moisture management.

Lined with Coolmax®

- Moves moisture away from your body
- Lining dries quickly
- Added performance keeps you cool and comfortable

Coolmax is a registered trademark of INVISTA (formerly DuPont Textile & Interiors) Airprene is a Trade Mark of Airprene LLC.

Thumb Spica - Neoprene



Features and Benefits

Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.

- Warmth and compression
- Removable flexible stay
- Restricts CMC and MCP joint
- Leaves IP joint free

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion
- Osteoarthritis
- Painful/unstable 1st mcp joint
- Rheumatoid arthritis
- Thumb strains
- Sprains and post-cast removal

*Elastic also available, please contact Customer Services.

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	14 - 15	TS/BE/S/L	TS/BE/S/R
Medium	15 - 16.5	TS/BE/M/L	TS/BE/M/R
Large	16.5 - 18	TS/BE/L/L	TS/BE/L/R
X-Large	18 - 19	TS/BE/XL/L	TS/BE/XL/R

Wrist and Thumb Spica - Neoprene



Features and Benefits

- Restricted wrist and thumb motion without limiting finger motion
- Easy fit and comfortable
- Lightweight and padded
- Contours to the anatomy of the wrist and thumb
- Breathable material
- Hook and loop fastening

Indications for Use

- Gamekeeper's thumb
- Post-cast healing
- Soft tissue injuries
- Selective scaphoid injuries

Material: Neoprene (breathable)

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	15 - 18	WTS/S/L	WTS/S/R
Medium	18 - 20	WTS/M/L	WTS/M/R
Large	20 - 23	WTS/L/L	WTS/L/R
X-Large	23 - 25	WTS/XL/L	WTS/XL/R

Wrist and Thumb Spica



Features and Benefits

- Two-piece dual compression with malleable aluminium thumb splint to hold thumb in abduction
- 5cm wide elasticated strap

Indications for Use

- Osteoarthritis
- Painful/unstable 1st mcp joint
- Rheumatoid arthritis

Material: Neoprene

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
X-Small	10 - 12.5	WTB/CX/BE/XS/L	WTB/CX/BE/XS/R
Small	12.5 - 15	WTB/CX/BE/S/L	WTB/CX/BE/S/R
Medium	15 - 17.5	WTB/CX/BE/M/L	WTB/CX/BE/M/R
Large	17.5 - 20	WTB/CX/BE/L/L	WTB/CX/BE/L/R
X-Large	20 - 22.5	WTB/CX/BE/XL/L	WTB/CX/BE/XL/R

Polypropylene Thumb Brace



Features and Benefits

Anatomically designed to provide comfortable and secure immobilisation of the thumb while preserving mobility of the other fingers and wrist.

Indications for Use

- DeQuervain's Syndrome
- Gamekeeper's thumb
- Tendonitis
- Thumb strains
- Sprains and post-cast removal

Manufactured from ultra-thin, rigid polyethylene with contact closure; includes two stockinette undersleeves.

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	15.2 - 17.8	37-5001	37-5011
Medium	17.8 - 20.3	37-5002	37-5012
Large	20.3 - 23	37-5003	37-5013
X-Large	23 - 25	37-5004	37-5014

CMC Brace



Features and Benefits

- Support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support.

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion

Material: Neoprene

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
X-Small	11.5 - 13.5	CMC/XSL	CMC/XSR
Small	14 - 15.5	CMC/SL	CMC/SR
Medium	16 - 19.5	CMC/ML	CMC/MR
Large	20 - 22	CMC/LL	CMC/LR
X-Large	22.5 - 24	CMC/XLL	CMC/XLR

1

2

3

4

5

6

7

8

9