

Diaped Duosoft Therapeutic



Features and Benefits

- Constructed with a multi-density material combination of plastazote and cellular urethane, these advanced insoles offer superb resilience, cushioning and shock absorption
- The insoles are designed to fit into both traditional footwear and 'additional depth' orthopaedic shoes
- The addition of an increased metatarsal support in the anatomical design is designed to facilitate off-loading of the metatarsal head area

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- Skin breakdown

UK Shoe Size	EU Shoe Size	Product Code
2 - 4.5	34 - 37	SB2335
5 - 6.5	36 - 38	SB2336
7 - 8.5	39 - 42	SB2337
9 - 10.5	43 - 45	SB2338
11 - 13	46 - 48	SB2339

Slim Tech



Features and Benefits

- The lightweight and easy fitting Slim Tech with its 3° low arch allows for heel contact and a good comfortable fit
- The orthotic is for use in ladies' dress shoes and sports boots
- Slip is reduced by the design of the hook heel

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code	UK Shoe Size	Product Code
3 - 4	SC3070	7 - 8	SC3072
5 - 6	SC3071	9 - 10	SC3073

Trio Slimfit



Features and Benefits

- The recessed medial heel post within the Trio Slimfit design offers control over pronation without risking irritation of the skin within the tighter fitting fashion shoe
- The forefoot contouring provides a stabilising effect upon the forefront joints and reduces cramping pain

Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density
X-Small	3 - 4.5	TSF-XS
Small	5 - 6.5	TSF-S
Medium	7 - 8.5	TSF-M
Large	9 - 10.5	TSF-L

Prostep



Features and Benefits

- Prostep is designed to restore the natural position of the talus by supporting the inside edge of the heel and supports the spring ligament
- Resists the valgus rotation of the rearfoot, which is associated with pronation
- Anti-bacterial top cover and forefoot and rearfoot integral posts of 0°, 3° & 5°
- Available in full or 3/4 length

Indications for Use

- Low arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	3/4 Medium	3/4 Firm	Full Medium	Full Firm
Kids Small	13.5 - 1.5	PSSLM-KS	PSSLF-KS	PSFLM-KS	PSFLF-KS
Kids Large	2 - 3	PSSLM-KL	PSSLF-KL	PSFLM-KL	PSFLF-KL
X-Small	3.5 - 4.5	PSSLM-XS	PSSLF-XS	PSFLM-XS	PSFLF-XS
Small	5 - 6	PSSLM-S	PSSLF-S	PSFLM-S	PSFLF-S
Medium	6.5 - 7.5	PSSLM-M	PSSLF-M	PSFLM-M	PSFLF-M
Large	8 - 9	PSSLM-L	PSSLF-L	PSFLM-L	PSFLF-L
X-Large	9.5 - 10.5	PSSLM-XL	PSSLF-XL	PSFLM-XL	PSFLF-XL
XX-Large	11 - 12	PSSLM-XXL	PSSLF-XXL	PSFLM-XXL	PSFLF-XXL

Trio



Features and Benefits

- The Trio insole is available in 3/4 and full length and is designed to resist excessive valgus rotation
- Inverts the calcaneus within the heel cup thereby maintaining mid tarsal joint
- Supports the proximal medial arch of the foot by providing navicular stability
- Each pair is supplied with additional 3° medial biplanar rearfoot posts and 3° lateral forefoot posts

Indications for Use

- Plantar fasciitis
- Heel pain
- Calf strain
- Anterior medial knee pain
- Hip instability

Size	UK Shoe Size	3/4 Length Medium Density	3/4 Length Firm Density	3/4 Length Extra Firm	Full Length Medium Density	Full Length Firm Density	Full Length Extra Firm
Kids Small	12.5 - 13.5	TSLM-KS	TSLF-KS	TSLE-KS	TFLM-KS	TFLF-KS	TFLE-KS
Kids Large	1 - 2.5	TSLM-KL	TSLF-KL	TSLE-KL	TFLM-KL	TFLF-KL	TFLE-KL
X-Small	3 - 4.5	TSLM-XS	TSLF-XS	TSLE-XS	TFLM-XS	TFLF-XS	TFLE-XS
Small	5 - 6.5	TSLM-S	TSLF-S	TSLE-S	TFLM-S	TFLF-S	TFLE-S
Medium	7 - 8.5	TSLM-M	TSLF-M	TSLE-M	TFLM-M	TFLF-M	TFLE-M
Large	9 - 10.5	TSLM-L	TSLF-L	TSLE-L	TFLM-L	TFLF-L	TFLE-L
X-Large	11 - 12.5	TSLM-XL	TSLF-XL	TSLE-XL	TFLM-XL	TFLF-XL	TFLE-XL
XX-Large	13 - 14.5	TSLM-XXL	TSLF-XXL	TSLE-XXL	TFLM-XXL	TFLF-XXL	TFLE-XXL

Control Tech Soft 3/4 Length



Features and Benefits

- The rearfoot wedging reduces the need for excessive arch height while maintaining support, comfort and providing superior shock absorption

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code	UK Shoe Size	Product Code
2 - 3	SC3055	7 - 8	SC3058
3 - 4	SC3056	9 - 10	SC3059
5 - 6	SC3057	11 - 12	SC3060

Equiflex



Features and Benefits

- By elevating the heel the Equiflex will restore ankle joint extension reducing the demand for compensation
- Lateral forefoot is supported at the level of the metatarsal heads to facilitate this and accommodate any additional postings that may be required
- Supplied with 0°, 3° and 5° rear posts and 6mm elevators

Indications for Use

- High arches
- Excessive pronation
- Forefoot pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density	Product Code Firm Density
X-Small	3.5 - 4.5	EQM-XS	EQF-XS
Small	5 - 6	EQM-S	EQF-S
Medium	6.5 - 7.5	EQM-M	EQF-M
Large	8 - 9	EQM-L	EQF-L
X-Large	9.5 - 10.5	EQM-XL	EQF-XL
XX-Large	11 - 12	EQM-XXL	EQF-XXL

Trio Easy Fit



Features and Benefits

- The Easy Fit is intrinsically a 3/4 length medially posted orthotic with an additional posting option for triplanar control of pronation featuring a lateral cutaway design
- This model maintains the essential features for pronation control, including an optional biplanar rearfoot post
- A medial forefoot post can be used to replace the first metatarsal cushion if required
- Comfort is achieved without compromising natural movement and comfort

Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density
Kids	1 - 2.5	TEF-K
X-Small	3 - 4.5	TEF-XS
Small	5 - 6.5	TEF-S
Medium	7 - 8.5	TEF-M
Large	9 - 10.5	TEF-L
X-Large	11 - 12.5	TEF-XL

Trio Heel Pain



Features and Benefits

- The Trio Heel Pain orthosis is contoured to support the navicular bone of the arch, directly reducing the cause of plantar fascia strain
- Higher loading forces under the heel can also lead to inflammation of the skin and supporting structures of the foot
- Features a large cushioning protection for the softer structure under the heel bones

Indications for Use

- Plantar fasciitis
- Plantar heel pain
- Calf strain
- Anterior medial knee pain
- Hip instability

Size	UK Shoe Size	Product Code Medium Density
Small	5 - 6.5	THP-S
Medium	7 - 8.5	THP-M
Large	9 - 10.5	THP-L

Control Tech Soft Full Length



Features and Benefits

- Available in low, moderate and high arch angles (4°, 6° & 8°)
- Plantar fascial grooves reduce arch irritation as well as provide comfort and superior shock absorption

Indications for Use

- Metatarsalgia
- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	4° Low Arch	6° Moderate Arch	8° High Arch
10 - 13	SC3086	-	-
1 - 2	SC3088	-	-
2 - 3	SC3090	SC3040	-
3 - 4	SC3091	SC3041	SC3031
5 - 6	SC3092	SC3042	SC3032
7 - 8	SC3093	SC3043	SC3033
9 - 10	SC3094	SC3044	SC3034
11 - 12	SC3095	SC3045	-
13 - 15	SC3096	SC3046	-

Interpod Diabetic



Features and Benefits

- Inverted rearfoot improves biomechanical function and enhances the windlass effect
- The 6° moderate arch height provides adequate functionality for the majority of patients
- PORON 96 top cover sets a new level of advanced underfoot technology
- The slow rebound cradles the foot and allows enhanced distribution of body weight over the foot's plantar surface
- Cradle functionality makes this the ultimate stock diabetic orthosis

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code	UK Shoe Size	Product Code
2 - 3	SC3130	9 - 10	SC3134
3 - 4	SC3131	11 - 12	SC3135
5 - 6	SC3132	13	SC3136
7 - 8	SC3133		

Control Tech Flex 3/4



Features and Benefits

- 3/4 length insole eliminates pronation moments about the subtalar joint created by the lateral heel cup
- Supplied with optional self adhesive covers
- Polypropylene integral posting over the subtalar joint axis

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	4° Low Arch	6° Moderate Arch	8° High Arch
3 - 4	SC3021	SC3011	SC3001
5 - 6	SC3022	SC3012	SC3002
7 - 8	SC3023	SC3013	SC3003
9 - 10	SC3024	SC3014	SC3004
11 - 12	SC3025	SC3015	SC3005

