


Universal VENUM Hip Brace



Features and Benefits

The one-size-fits-all VENUM Soft Hip Abduction brace incorporates the most advanced technology and materials. A revolutionary soft hip management system for patients of all different shapes and sizes providing maximum comfort and support. The LSO part of the brace operates utilising a cincher system to maximise control. The belt thigh cuff and hinge adjust to fit the dimension of most patients. The brace applies a dynamic abduction force, so that the femoral head remains in the socket.

- Universal design
- Cincher technology
- Size indicator displays belt size to ensure correct and comfortable fitting
- Joint with adjustable abduction/adduction, flexion/extension
- Single or double-pull cable pulley system
- Stock reduction, saves money and space
- Materials are comfortable, breathable, soft, wick-away perspiration and are durable and resilient
- Waist circumference range: 60cm - 157.5cm
- Extender panel adds an additional 25cm (sold separately)
- Ergonomic arthritic grips can be customised to patient size

Indications for Use

- For patients at risk of hip dislocation, suitable for all levels of instability

Material: 30% Polyurethane 70% Polyamide

Size	Description	Product Code
Universal	VENUM Soft Hip Abduction Brace	VHHIPABD
Universal	Long Joint*	TLCHIPJNTLG


* Note: The Venum Soft Hip Aduction Brace comes with a standard joint for patients of height below 188cm, for patients over 188cm a long joint will be required.

VENUM Universal Hip Brace Size Chart

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5		

To fit the VENUM Hip Brace correctly, please choose the desired size on the slide panel from the chart.

Hip Control System



Features and Benefits

- Preformed modular plastic waist and thigh components
- Combination Kinesio joint, controls flexion/extension, abduction/adduction
- Hip and waist components supplied with comfort liners
- Component sizes can be mixed for optimum fit
- Kinesio offset hinge with abduction/adduction settings (7.5° increments) and flexion/extension settings (from 0° to 120° in 15° increments)

Indications for Use

- This modular hip orthosis system offers functional and effective support in the treatment of various hip pathologies

Hip / Waist Component					Kinesio Hinge		
Size	Waist Circ. (cm)	Hip Circ. (cm)	Product Code Left	Product Code Right	Size	Height (cm)	Product Code
Small	60 - 82	68 - 90	M.2122	M.2112	Regular	< 180	M.2303
Medium	94 - 115	90 - 110	M.2123	M.2113	Long	> 181	M.2304
Large	39 - 44	106 - 126	M.2124	M.2114			
X-Large	99 - 120	111 - 131	M.2125	M.2115			

Thigh Component				
Size	Proximal Circ. (cm)	Distal Circ. (cm)	Product Code Left	Product Code Right
Small	44 - 53	34 - 41	M.2222	M.2212
Medium	48 - 60	36 - 45	M.2223	M.2213
Large	57 - 68	42 - 49	M.2224	M.2214
X-Large	62 - 76	46 - 56	M.2225	M.2215