

VENUM Universal TLSO PPX



Features and Benefits

The VENUM TLSO PPX is the ultimate one-size-fits-all semi rigid comfort brace for the thoraco lumbar-sacral region. Primarily, this brace reduces pain, whilst providing support and immobilisation to the thoraco lumbar-sacral spine. The brace features unique anterior pads and in combination with the highly effective cincher system is designed to produce maximal intra-abdominal pressure to distract the vertebral bodies of the spine. As the brace is universal in nature, it benefits from fitting every patient type and can be stored easily where limited storage is available.

- Reduced stock holding
- Cincher technology
- Cost effective
- Step-down, streamline design
- Motion restriction of the thoracic and lumbar regions L5 - T4
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic arthritic grips can be customised to patient size by extending or retracting cable length
- Made from durable, breathable, soft, hypo-allergenic materials that wicks away perspiration
- Fits waist circumferences from 60cm - 157.5cm; optional extender panel adds 25cm to circumference

Indications for Use

- Chronic back pain
- Degenerative, herniated or bulging disc
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- Spinal instability
- Burst Fx
- Pre or post-op spinal thoraco-lumbar stabilisation
- Pre or post-op disk surgery

Size	Description	Product Code
Universal	VENUM TLSO PPX	VIPHYBTLSO4PPX
Universal	Extender Panel	STEXT

VENUM Universal HYBRID TLSO



Features and Benefits

The VENUM HYBRID TLSO is the ultimate one-size-fits-all semi rigid comfort brace for the thoraco lumbar-sacral region. This brace has the ability to contour to a kyphotic-shaped spine and is more forgiving in the treatment of elderly patients who would struggle with the rigidity of a TLSO with anterior extension. Primarily, this brace reduces pain, whilst providing support and immobilisation to the thoraco lumbar-sacral spine. The brace features unique anterior pads and in combination with highly effective cincher system is designed to produce maximal intra-abdominal pressure to distract the vertebral bodies of the spine. As the brace is universal in nature, it benefits from fitting every patient type and can be stored easily where limited storage is available.

- Reduced stock holding
- Cincher technology
- Cost effective
- Step-down design
- Streamline design
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic arthritic grips can be customised to patient size by extending or retracting cable length
- Made from durable, breathable, soft, hypo-allergenic materials that wick away perspiration
- Fits waist circumferences from 60cm - 157.5cm; optional extender panel adds 25cm to circumference

Indications for Use

- Chronic back pain
- Kyphosis
- Osteoporosis or osteoporotic fractures
- Degenerative, herniated or bulging disc
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- Spinal instability
- Burst Fx
- Pre or post-op spinal lumbar stabilisation
- Pre or post-op disk surgery

Size	Description	Product Code
Universal	VENUM HYBRID TLSO	VIPHYBTLSO3DLK
Universal	Extender Panel	STEXT