

Fracture Walker Brace



Features and Benefits

- Impact-resistant nylon foot piece with removable foam liners
- Liners attach to uprights with touch hook fastening
- Liners are secured to uprights with broad circumferential touch and close fastening
- Universal left or right

Indications for Use

- Post-op ankle stabilisation
- Ankle fracture
- Stress fractures
- Soft tissue and ligament injuries

Size	UK Shoe Size	Product Code
Small	Up to 5	FW17/S
Medium	5 - 8	FW17/M
Large	8 - 11	FW17/L
X-Large	12 +	FW17/XL

Rebound Air Walker



Features and Benefits

The Rebound Air Walker features linked medial and lateral air bladders with a single inflation point to provide balanced compression. These air bladders are shaped to cradle the achilles, holding the heel in place for extra stability.

- Moisture wicking soft goods and ventilated panels prevent microbial build up and helps to keep the skin dry and comfortable
- Straps fasten in front for easy application, adjustment and removal
- EVA material provides increased shock absorption for the heel
- Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression
- Flex Edge over moulding conforms to each patient's unique anatomy to help limit the formation of pressure points
- Rocker bottom promotes a natural, stable gait, whilst a semi-soft toe bumper that can be trimmed as required provides added protection

Indications for Use

- Grade 2 and 3 ankle sprains
- Stable fractures
- Post-operative use
- Trauma
- Rehabilitation use

Size	Men's Shoe Size	Ladies' Shoe Size	Product Code
X-Small	2 - 4	3.5 - 5.5	B-242900001E
Small	4 - 7	5.5 - 8.5	B-242900002E
Medium	7 - 10.5	8.5 - 12	B-242900003E
Large	10.5 - 12.5	12 - 14	B-242900004E
X-Large	12.5 - 15	14 - 16.5	B-242900005E

Rebound Diabetic Walker



Features and Benefits

Uniquely designed for effective treatment of plantar foot ulcers, the Rebound Diabetic Walker enables easy access and re-fitting during wound care, whilst improving the patient's mobility during the healing process by including a rocker bottom outsole for smooth roll-over.

- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote top layer that evenly distributes the pressure
- Anatomically designed shell that minimises shear forces by immobilising the foot and ankle and 90 degrees whilst providing protection
- Flex Edge technology for relief of pressure points, along with a patent-pending pressure relief valve to prevent over-inflation
- Positioning portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Soft and breathable machine washable liner for hygiene promotion and prevention of infection

Indications for Use

- Plantar foot ulcers

Size	Men's Shoe Size	Ladies' Shoe Size	Product Code
Small	4.5 - 7	5.5 - 8.5	B-242500002
Medium	7 - 10.5	8.5 - 12	B-242500003
Large	10.5 - 12.5	12 - 14	B-242500004
X-Large	12.5 - 15	14 - 16.5	B-242500005

VACOped®



Features and Benefits

Inspired by ski boots, the VACOped features a unique design to enable patients to weight-bear earlier, enabling a faster return to daily activities.

- Inner vacuum cushion provides equal and firm support by adjusting itself precisely to the shape of every foot and lower leg and allowing adjustments to be made for swelling
- Adjustable ROM bar allows the boot to be fixed or permit movement dependent on the indication
- Durable, yet lightweight German-engineered plastic shell provides the injury site with the necessary stability and allows an increased freedom of movement in comparison to a plaster cast
- Multiple washable and exchangeable liners promote hygiene

Indications for Use

- Achilles tendon ruptures
- Ankle fractures
- Metatarsal fractures
- Calcaneus fractures
- Ankle replacements
- Arthrodesis

Size	UK Shoe Size	Product Code
Small	3 - 6.5	VPED/S
Medium	6.5 - 9	VPED/M
Large	9 - 10+	VPED/L

Repose Foot Protector



Features and Benefits

Repose Foot Protectors are a unique product designed specifically to minimise the risk of pressure damage to heels. A significant proportion of pressure ulcers occur on the heels due to a combination of pressure and friction/shear. Whilst support surfaces assist in reducing pressure, only a specific device such as the Repose Foot Protector can ensure that pressure on the foot is relieved totally. Clinical studies have shown these heel protectors to be effective in a wide variety of settings.

- May be used in bed or with feet elevated, easily secured with stocking or light bandage
- Can be used with other dynamic or static support surfaces
- Reduces impact of pressure, friction and shear
- Maximum patient weight = 139kg

Indications for Use

- Pressure relief

Size	Length (cm)	Width (cm)	Height (cm)	Product Code
Small	38	20	18	6501100
Large	51	28	18	6021100

Repose Wedge



Features and Benefits

The Repose Wedge has been designed primarily for off-loading pressure from the heels and can be used wherever support or elevation of the lower legs is required. The Repose Wedge offers the same pressure reducing properties as other Repose products. Particularly suited to leg elevation to aid vascular conditions and a useful aid in posture support.

Indications for Use

- The Repose Wedge is a versatile product suitable for a wide range of applications within pressure area

Size	Product Code
Universal (elastic strap)	6711111
Universal (clear strap)	6700000