

Comfy Air Hand Orthosis



Features and Benefits

- The Comfy Air Hand Orthosis is an exceptionally useful splint for the severely contracted hand
- It has an adjustable air bladder that allows gradual extension of the MP, DIP and PIP joints, as well as an adjustable wrist for customised fit
- Ingeniously designed by a therapist, one simply inserts the un-fitted splint with the bladder under the fingers and then pumps air to inflate the bladder and increase finger extension. (Headliner material only)
- Ambidextrous

Indications for Use

- Ideal for the severely contracted hand
- Inflatable air bladder allows for gradual MP, DIP and PIP joint extension
- Adjustable wrist
- Provides low pressure stretch therapy
- Hand air bladder can be used separately without full orthosis as needed
- Can fit into even the most tightly clenched hand

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	HA101-A

Short Wrist Brace



Features and Benefits

Standard single compression one-piece wrist brace provides good conformity around the wrist and hand.

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

Material: Neoprene

		Product Code - Neoprene Beige	
Size	Wrist Circ. (cm)	Left	Right
X-Small	10 - 12.5	SWB/E/XS/L	SWB/E/XS/R
Small	12.5 - 15	SWB/E/S/L	SWB/E/S/R
Medium	15 - 17.5	SWB/E/M/L	SWB/E/M/R
Large	17.5 - 20	SWB/E/L/L	SWB/E/L/R
X-Large	20 - 22.5	SWB/E/XL/L	SWB/E/XL/R

Short Wrist Brace - Elastic



Features and Benefits

Standard single compression one-piece wrist brace provides good conformity around wrist and hand.

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

Material: Elastic

		Product Code - Elastic Beige	
Size	Wrist Circ. (cm)	Left	Right
X-Small	10 - 12.5	SWB/XSL	SWB/XSR
Small	12.5 - 15	SWB/SL	SWB/SR
Medium	15 - 17.5	SWB/ML	SWB/MR
Large	17.5 - 20	SWB/LL	SWB/LR
X-Large	20 - 22.5	SWB/XLL	SWB/XLR

Long Wrist Brace



Features and Benefits


- Longer two-piece design giving support and comfort for continuous use
- Shape allows increased finger mobility and minimal thumb web discomfort
- Fully bound edges

Indications for Use

- Flaccid drop wrist
- Rheumatoid arthritis
- Ulnar and radial stability

Material: Coolmax®

Size	Wrist Circ. (cm)	Product Code - Coolmax® Beige	
		Left	Right
X-Small	Up to 14	LWB/CX/BE/XS/L	LWB/CX/BE/XS/R
Small	14 - 16	LWB/CX/BE/S/L	LWB/CX/BE/S/R
Medium	16.5 - 19	LWB/CX/BE/M/L	LWB/CX/BE/M/R
Large	19 - 21.5	LWB/CX/BE/L/L	LWB/CX/BE/L/R
X-Large	Over 21.5	LWB/CX/BE/XL/L	LWB/CX/BE/XL/R



About Coolmax®

CoolMax is the industry leader when it comes to moisture management.

Lined with Coolmax®

- Moves moisture away from your body
- Lining dries quickly
- Added performance keeps you cool and comfortable

Coolmax is a registered trademark of INVISTA (formerly DuPont Textile & Interiors) Airprene is a Trade Mark of Airprene LLC.

Thumb Spica - Neoprene



Features and Benefits

Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free.

- Warmth and compression
- Removable flexible stay
- Restricts CMC and MCP joint
- Leaves IP joint free

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion
- Osteoarthritis
- Painful/unstable 1st mcp joint
- Rheumatoid arthritis
- Thumb strains
- Sprains and post-cast removal

*Elastic also available, please contact Customer Services.

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	14 - 15	TS/BE/S/L	TS/BE/S/R
Medium	15 - 16.5	TS/BE/M/L	TS/BE/M/R
Large	16.5 - 18	TS/BE/L/L	TS/BE/L/R
X-Large	18 - 19	TS/BE/XL/L	TS/BE/XL/R