

Wrist and Thumb Spica - Neoprene



Features and Benefits

- Restricted wrist and thumb motion without limiting finger motion
- Easy fit and comfortable
- Lightweight and padded
- Contours to the anatomy of the wrist and thumb
- Breathable material
- Hook and loop fastening

Indications for Use

- Gamekeeper's thumb
- Post-cast healing
- Soft tissue injuries
- Selective scaphoid injuries

Material: Neoprene (breathable)

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	15 - 18	WTS/S/L	WTS/S/R
Medium	18 - 20	WTS/M/L	WTS/M/R
Large	20 - 23	WTS/L/L	WTS/L/R
X-Large	23 - 25	WTS/XL/L	WTS/XL/R

Wrist and Thumb Spica



Features and Benefits

- Two-piece dual compression with malleable aluminium thumb splint to hold thumb in abduction
- 5cm wide elasticated strap

Indications for Use

- Osteoarthritis
- Painful/unstable 1st mcp joint
- Rheumatoid arthritis

Material: Neoprene

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
X-Small	10 - 12.5	WTB/CX/BE/XS/L	WTB/CX/BE/XS/R
Small	12.5 - 15	WTB/CX/BE/S/L	WTB/CX/BE/S/R
Medium	15 - 17.5	WTB/CX/BE/M/L	WTB/CX/BE/M/R
Large	17.5 - 20	WTB/CX/BE/L/L	WTB/CX/BE/L/R
X-Large	20 - 22.5	WTB/CX/BE/XL/L	WTB/CX/BE/XL/R

Polypropylene Thumb Brace



Features and Benefits

Anatomically designed to provide comfortable and secure immobilisation of the thumb while preserving mobility of the other fingers and wrist.

Indications for Use

- DeQuervain's Syndrome
- Gamekeeper's thumb
- Tendonitis
- Thumb strains
- Sprains and post-cast removal

Manufactured from ultra-thin, rigid polyethylene with contact closure; includes two stockinette undersleeves.

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	15.2 - 17.8	37-5001	37-5011
Medium	17.8 - 20.3	37-5002	37-5012
Large	20.3 - 23	37-5003	37-5013
X-Large	23 - 25	37-5004	37-5014