


Oasis TLSO



Model wearing PPK version shown above.



Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid LSO braces.

Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates.

The Oasis TLSO and TLSO LP 3 Panel DLK includes a Dorsal Lumbar Kit (DLK) with Sternal Panel to prevent forward flexion. Effective within T4 - S1.

Indications for use:

- Burst Fx
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Kyphosis
- Laminectomy
- Osteoporosis
- Paraplegia
- Post-Operative
- Spinal Stabilisation
- Spinal Stenosis
- Spondylolysis
- Spondylolisthesis

Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OATL4-XS-PPK
Small	76.2 - 91.4	88.9 - 104.14	OATL4-S-PPK
Medium	86.4 - 101.6	99 - 114.3	OATL4-M-PPK
Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-PPK
X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-PPK
XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-PPK
XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-PPK

Available in SPK & PPK versions, other versions available to order.

Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OATL4-XS-SPX
Small	76.2 - 91.4	88.9 - 104.14	OATL4-S-SPX
Medium	86.4 - 101.6	99 - 114.3	OATL4-M-SPX
Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-SPX
X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-SPX
XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-SPX
XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-SPX

Steeper Group
Unit 3, Stourton Link
Intermezzo Drive
Stourton Link
Leeds
LS10 1DF

Tel: +44 (0) 870 240 4133
Email: customerservices@steepergroup.com

www.steepergroup.com