

We hope you are happy with your adapted or repaired footwear.

Our goal at Steeper is to always achieve the best possible outcomes from the orthotic devices we manufacture, supply or adapt for patients referred into our care.

# Footwear Repairs & Adaptations Care & Usage Instructions

The footwear is intended for the single exclusive use by the named patient.

## Repairs & adaptations to footwear

1. Please take care when removing footwear that has been adapted.
2. Keep all footwear away from direct sources of heat.
3. In the event that you experience any problems to your adapted footwear or its construction, please contact the hospital department where they were issued.
4. Your adapted footwear should be periodically checked for wear and reported to your orthotist.
5. An annual review of your adapted footwear is advisable by the orthotist at your hospital.

Ref No: .....

Orthotist: .....

Patient Name: .....

Unit 3, Stourton Link  
Intermezzo Drive  
Leeds, LS10 1DF  
Tel: +44 (0) 113 207 0435  
Email: [orthoticenquiries@steepergroup.com](mailto:orthoticenquiries@steepergroup.com)

[www.steepergroup.com](http://www.steepergroup.com)

