

We trust you are satisfied with your new insoles and they prove to be effective.

Our goal at Steeper is to always achieve the best possible outcomes from the orthotic devices we manufacture, supply or adapt for patients referred into our care.

Insoles Care & Usage Instructions

This orthosis is intended for the single exclusive use by the named patient.

Custom-made device

1. Always use your insoles in correctly fitting sensible footwear. Your orthotist can advise on appropriate footwear.
2. Wear insoles for short periods initially.
3. Clean insoles by wiping with a damp cloth.
4. If you experience any problems with the fit of your insoles, please contact the hospital department where they were issued.
5. Your Orthotist will advise you on the correct review process

Ref No:

Orthotist:

Patient Name:

Steeper Group, Unit 3, Stourton Link,
Intermezzo Drive, Leeds, LS10 1DF
Tel: 0113 207 0435
Email: orthoticenquiries@steepergroup.com

www.steepergroup.com

