

Wear & Care Instructions



Padded PF Splint

Steeper manufactures a quality range of orthotic products and to help you care for the product, this leaflet provides information on how and when to use your prescribed orthosis, how to care for it and how to contact us if you have a problem or require a replacement.

Care & Cleaning

Close all contact closures and hand wash in warm water with mild soap. Product should be airdried, do not use high temperatures to wash or dry. Remove the fabric cover to wash.

Applying the Orthosis

1. Place affected foot inside the boot ensuring the heel is in position.
2. Wrap the straps over the product onto opposite side and attach straps onto hook portion of the product. Straps should be pulled gently for comfortable fit.
3. Straps should be pulled gently for comfortable fit.

IMPORTANT

The advice and information contained in this document is not intended as a substitute for medical counselling. If you have any questions about treatment for a specific condition, please check with a health care professional.



Unit 3 Stourton Link,
Intermezzo Drive,
Leeds LS10 1DF
Tel: +44 (0) 113 207 0435
Email: orthoticenquiries@steepergroup.com

www.steepergroup.com

 **steeper**
Creating life's turning points, together