

PeroSupport.tec

Indications for Use

- Peroneal paralysis, weakness of ankle dorsiflexors, inflammation of the ankle joint

Features and Benefits

- Unique, single-handed magnetic fastening system - suitable for hemiparesis
- Great energy storage and return
- Limits the plantar flexion in the upper ankle joint and therefore causes the foot to lift during swing phase
- Supplied with two sets of cushioning pads
- Can be trimmed with sturdy scissors to size

Size	UK Shoe Size	Height	Product Code	
			Left	Right
X-Small	2-3.5	32.5cm	197R4/L0	197R4/R0
Small	3-6	33.5cm	197R4/L1	197R4/R1
Medium	6-8	34.5cm	197R4/L2	197R4/R2
Large	8-10	35.5cm	197R4/L3	197R4/R3
X-Large	10-13	37.5cm	197R4/L4	197R4/R4



Unique fastening system

PeroSupport.pro

Indications for Use

- Peroneal paralysis, weakness of ankle dorsiflexors, inflammation of the ankle joint

Features and Benefits

- Great energy storage and return
- Limits the plantar flexion
- Prevents lateral ankle twisting
- Prevents side-slips
- Supplied with two sets of cushioning pads
- Can be trimmed with sturdy scissors to size

Size	Shoe Size	Height	Product Code	
			Left	Right
Small	3-6	34cm	197R6/L1	197R6/R1
Medium	6-8	35cm	197R6/L2	197R6/R2
Large	8-10	35cm	197R6/L3	197R6/R3
X-Large	10-13	36cm	197R6/L4	197R6/R4

