


## Structural and Socket Components - 30mm Heavy Duty Titanium and Stainless Steel




**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

30mm Heavy Duty Titanium Structural Components			
Part Number	Manufacturer	Description	Weight Limit
1 NAAASMP450	TL	Socket Mounting Plate - includes expulsion valve insert and fabrication components. This lamination plate is supplied with an expulsion valve as standard, but can be used without, if required.	205kg
2 NMA10A3/T	SF	Heavy Duty Socket Adapter	175kg
3 NMA10A90/T	SF	Heavy Duty Socket Adapter	175kg
4 NMA11A11/T	SF	Heavy Duty Tube Clamp Adapter	175kg
5 NMA15A4/T225	SF	Heavy Duty 214mm Tube Adapter	175kg
NMA15A4/T445	SF	Heavy Duty 414mm Tube Adapter	175kg
6 NMA10A42/T	SF	Heavy Duty Slide Adapter	175kg

**Note:** The weight limit is 168kg for the heavy duty aluminium components above for patients of K4 activity level.

Manufacturer: TL= Trulife SF = Streifeneder



Heavy Duty Stainless Steel College Park Structural Components			
Part Number	Manufacturer	Description	Weight Limit
1 NHD-TTK-PA	CP	Transtibial Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 250mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter)	181kg
2 NHD-TFK-PA	CP	Transfemoral Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 430mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter)	181kg
3 NLP-RHD-S	CP	4-Hole Lamination Plate, Round, Heavy Duty, SS	181kg

Manufacturer: CP = College Park

## Structural and Socket Components - 30mm & 34mm Heavy Duty Titanium



**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

30mm & 34mm Heavy Duty Titanium Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NAAASMP450	TL	Socket Mounting Plate - includes expulsion valve insert and fabrication components. This lamination plate is supplied with an expulsion valve as standard, but can be used without, if required.	205kg
2	NHD-425BH	APC	Heavy Duty Titanium 34mm Pyramid Adapter	227kg
3	NHD-AL-PYLON	APC	34mm Aluminum Heavy Duty Pylon	227kg

**Note:** The weight limit is 168kg for the heavy duty aluminium components above for patients of K4 activity level.

Manufacturer: TL = Trulife APC = American Prosthetic Components



34mm Heavy Duty Titanium Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NHD-100C	APC	Heavy Duty Titanium Tube Clamp	227kg
2	NHD-110LONG	APC	Heavy Duty Titanium Receiver/Pylon Adapter - Long (82mm) Pylon	227kg
3	NHD-105C	APC	Heavy Duty Titanium Tube Clamp with 5mm Offset	227kg
4	NHD-132D	APC	Heavy Duty Titanium Dual Receiver Clamp	227kg
5	NHD-740L	APC	Heavy Duty Titanium 4-Prong Lamination Pyramid	227kg
6	NHD-100C-A	APC	Heavy Duty Tube Clamp with 30mm Receiver	227kg

Manufacturer: APC = American Prosthetic Components