

## Structural and Socket Components - 30mm & 34mm Heavy Duty Titanium

### 30mm Trulife Components Weight Limit - 205kg (unless otherwise stated)

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.

#### 30mm Heavy Duty Titanium Structural Components

Part Number	Description	Weight Limit
1 NAAAHTI210	4-Hole Male Pyramid Adapter	205kg
2 NAAAHTI217	Splined 4-Hole Male Rotatable Pyramid Adapter	205kg
3 NAAAHTI218	Splined 4-Hole Female Rotatable Receiver Adapter	205kg
4 NAAAHTI221	35mm Long Double Ended Adapter	204kg
5 NAAAHTI225	Clamp Adapter	204kg
6 NAAAHTI200	Bonded Pylon 420mm	205kg
NAAASMP350	Socket Mounting Plate - Includes Expulsion Valve Insert and fabrication components. This lamination plate is supplied with an expulsion valve as standard, but can be used without, if required.	160kg
7 NAAASMP450	Socket Mounting Plate - Includes Expulsion Valve Insert and fabrication components. This lamination plate is supplied with an expulsion valve as standard, but can be used without, if required.	205kg

Note: The weight limit is 168kg for the heavy duty aluminium components above for patients of K4 activity level.

### 34mm APC Components Weight Limit - 227kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.

#### 34mm Heavy Duty Titanium Structural Components

Part Number	Description	Weight Limit
1 NHD-100C	Heavy Duty Titanium Tube Clamp	227kg
2 NHD-110LONG	Heavy Duty Titanium Receiver / Pylon Adapter - Long (82mm) Pylon	227kg
3 NHD-105C	Heavy Duty Titanium Tube Clamp with 5mm Offset	227kg
4 NHD-132D	Heavy Duty Titanium Dual Receiver Clamp	227kg
5 NHD-740L	Heavy Duty Titanium 4-Prong Lamination Pyramid	227kg

