


Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Titanium 30mm Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NTI-100C-FL	APC	Titanium Tube Clamp Adapter	166kg
2	NTI-100C-SH	APC	Titanium 30mm Tube Clamp Adapter Short - This adapter has been designed to provide a build that is approximately 10mm shorter than the standard clamp, with the knurled inner surface of the clamp providing the necessary grip to prevent rotation	166kg
3	NTI-205C	APC	Titanium 5 Degree Angled Tube Clamp Adapter (pictured)	166kg
	NTI-210C		Titanium 10 Degree Angled Tube Clamp Adapter	166kg
	NTI-215C		Titanium 15 Degree Angled Tube Clamp Adapter	166kg
	NTI-220C		Titanium 20 Degree Angled Tube Clamp Adapter	166kg
4	NTI-120P	APC	Titanium Tube Clamp w/ Pyramid Adapter	166Kg
5	NTI-740L	APC	Titanium 4 Prong Laminating Base w/ Fixed Pyramid Adapter	166Kg
6	NTI-110LONG	APC	Titanium Pylon with Receiver Adapter (Overall length = 96.5mm) - This pylon, with its integral female receiver has an overall length of nearly 100mm, but can be cut down very short and with the use of NTI-100C-SH can provide a double adapter of less than 60mm, with rotation.	166Kg
	NTI-110SH	APC	Alternatively, use the NTI-110SH (pictured) which is approximately 50mm in length.	
7	NTI-105C	APC	Titanium Heavy Duty Tube Clamp, 5mm Offset	166kg
8	NTI-125P	APC	Titanium Pyramid / Receiver Adapter	166kg
9	NTI-112-P	APC	Titanium 30mm Pylon (SHORT) w/Pyramid Adapter - Overall Length of 2.46" Accommodates 1 - NTI-100C	166kg
10	NTI-130D	APC	Titanium 30mm Dual Receiver Adapter	166kg
	NTI-132D	APC	Titanium 32mm Dual Receiver Adapter	166kg
	NTI-138D	APC	Titanium 38mm Dual Receiver Adapter	166kg
	NTI-145D	APC	Titanium 45mm Dual Receiver Adapter	166kg
	NTI-152D	APC	Titanium 52mm Dual Receiver Adapter	166kg
	NTI-160D	APC	Titanium 60mm Dual Receiver Adapter	166kg
	NTI-168D	APC	Titanium 68mm Dual Receiver Adapter	166kg
	NTI-175D	APC	Titanium 75mm Dual Receiver Adapter	166kg

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Note: See pages 246-247 for the Mix and Match APC Alignment Adapters Component Build Height Information

Manufacturer: APC = [American Prosthetic Components](#)



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Titanium 30mm Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NTI-465TRF	APC	Titanium 4 Hole Clamping Adapter (Centre Thread) with Rotatable Receiver	166kg
2	NTI-460TRM	APC	Titanium 4 Hole Clamping Adapter (Centre Thread) with Rotatable Pyramid	166kg
3	NTI-460TB	APC	Titanium 4 Hole Clamping Adapter with Centre Thread - This part number is only for the baseplate used on NTI-465TRF and NTI-460TRM. It will obviously take either NTI-50FS or NTI-00MS or any component with an M36 threaded top.	166kg
4	NTI-425BH	APC	Titanium 4 Hole Pyramid Adapter with Centre Hole	166kg
5	NTI-400B	APC	Titanium 4 Hole Pyramid Adapter	166kg
6	NTI-400A	APC	Titanium Offset 4 Hole Pyramid Adapter (Axial) - This adapter has the pyramid offset by 5mm in one direction.	166kg
7	NTI-400D	APC	Titanium Offset 4 Hole Pyramid Adapter (Diagonal) - This adapter has the pyramid offset by 5mm in two directions.	166kg
8	NTI-00MS	APC	Titanium Rotatable Pyramid with Male Thread and Centre Hole - This item can be used in the baseplate NTI-460TB, NSS-800BASE or and adapter with an M36 thread.	166kg
9	NTI-50FS	APC	Titanium Receiver with Male Thread - This item can be used in the baseplate NTI-460TB or any adapter with an M36 thread.	166kg
10	NTI-60FS	APC	Titanium Rotatable Receiver with Female Thread Height 36.3mm	166kg
11	NTI-01MS	APC	Titanium Rotatable Pyramid with Threaded Clamp	166kg

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Note: See pages 246-247 for the Mix and Match APC Alignment Adapters Component Build Height Information

Manufacturer: APC = [American Prosthetic Components](#)

Structural and Socket Components - 30mm Titanium & Stainless Steel

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.



Titanium 30mm Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NTI-17TC2-BLACK	APC	440mm Pylon with Titanium Receiver Adapter BLACK	166kg
2	NTI-10TC2-BLACK	APC	262mm Pylon with Titanium Receiver Adapter BLACK	166kg
3	NTI-10TC2-I-BLACK	APC	262mm Pylon with Internal Titanium Receiver Adapter BLACK	166kg

This pylon has its adapter bonded to the inside of the tube, which some clinicians prefer the appearance of and also allows it to be cut much shorter, should the need arise.

Note: Whilst the components showing the blue water symbol are therefore waterproof, these components do not have drainage slots.



Stainless Steel 30mm Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NSS-800LM	APC	SS 3 Prong Adapter with Rotatable Pyramid	125kg
2	NSS-850LM	APC	SS 3 Prong Adapter with Rotatable Receiver	125kg
	NSS-800BASE	APC	SS 3 Prong Adapter Threaded with Female Thread This part number is only for the adapter, as used on NSS-800LM and on NSS-850LM with their stainless steel pyramid and female receiver. It will take either NTI-50FS or NTI-00MS to provide a lighter build, or any component with an M36 threaded top.	125kg
3	NSS-50FS	APC	Stainless Steel Receiver with Male Thread This item can be used in the baseplate NTI-460TB or any adapter with an M36 thread.	136kg
4	NTI-420L	APC	Titanium Round 4-Hole Laminating Plate This simple 4.8mm thick lamination plate, with 4 M6 threaded holes, comes with 4 of each M6 x 20mm and M6 x 25mm water resistant screws and is ideal for laminate sockets, but can also be used on thermoplastic sockets, provided a flat surface can be achieved onto which a 4 hole adapter can be bolted.	166kg
5	NSS-01MS	APC	Stainless Steel Rotatable Pyramid with Threaded Clamp	136Kg
6	NSS-00MS	APC	Rotatable Pyramid with Male Thread	136kg
7	NTI-TF900SET	APC	Titanium Thermoplastic Socket Attachment Plate Set This two-piece socket attachment plate has been designed for use on thermoplastic sockets, to provide a flat mounting surface for any 4-hole adapter, but can be used on laminated sockets as well.	166Kg
8	NAL-TF900SET	APC	APC Aluminium Thermoplastic Socket Attachment Plate Set This two piece socket attachment plate has been designed for use on thermoplastic sockets, to provide a flat mounting surface for any 4 hole adapter, but can be used on laminated sockets as well.	100Kg
9	NAL-420L	APC	APC Aluminium Round 4-Hole Laminating Plate Threaded This simple aluminium plate, with 4 threaded holes, is ideal for use with laminate sockets, but can also be used on thermoplastic sockets, provided a flat surface can be achieved onto which a 4 hole adapter can be bolted.	100Kg

Note: These Stainless Steel components are not waterproof or water resistant.
Note: See pages 246-247 for the Mix and Match APC Alignment Adapters Component Build Height Information
Manufacturer: APC = American Prosthetic Components