



**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

| 34mm Heavy Duty Stainless Steel College Park Structural Components |              |             |  |       |
|--|--------------|-------------|--|-------|
| Part Number  | Manufacturer | Description | Weight Limit   |       |
| 1  | N4HPA-S      | CP          | 4-Hole Pyramid Adapter, Standard   | 181kg |
| 2  | N4HPR-S      | CP          | 4-Hole Pyramid Receiver, Standard  | 181kg |
| 3  | NTA-S250-34  | CP          | Low-Profile Tube Adapter, 250mm length, 34mm diameter  | 181kg |
| 4  | NTA-S430-34  | CP          | Low-Profile Tube Adapter, 430mm length, 34mm diameter  | 181kg |
| 5  | NTC-S-34     | CP          | Tube Clamp, 34mm diameter  | 181kg |
| 6  | NHD-TTK-PA   | CP          | Transtibial Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 250mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter)  | 181kg |
| 7  | NHD-TFK-PA   | CP          | Transfemoral Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 430mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter) | 181kg |

Manufacturer: CP = College Park