


Foam Cosmesis



Transfemoral Cosmetic Soft Foam Cover
Large circumferences with 2-step-bore (Ø 60mm/Ø 32mm)

Part Number	Side	Calf Circ.	Knee Height
1 NMA65A10	Left/Right	52cm	52cm

Cosmetic Soft Foam Cover
Modular - transfemoral prosthesis with 2-step-bore (Ø 60mm/ Ø 32mm)

Part Number	Side	Calf Circ.	Knee Height
2 NMA65A11/44L	Left	44cm	52cm
NMA65A11/44R	Right	44cm	52cm

Cosmetic Soft Foam Cover
Modular - transfemoral prosthesis with conical 3-step-bore (Ø 70mm/Ø 50mm/Ø 32mm)

Part Number	Side	Calf Circ.
3 NMA65A12/40L	Left	40cm
NMA65A12/40R	Right	40cm
NMA65A12/44L	Left	44cm
NMA65A12/44R	Right	44cm

Cosmetic Soft Foam Cover
Modular - transfemoral prosthesis with conical 3-step-bore (Ø 60mm/Ø 50mm/Ø 32mm)

Part Number	Side	Calf Circ.	Knee Height
4 NMA65A13/40L	Left	40cm	52cm
NMA65A13/40R	Right	40cm	52cm
NMA65A13/44L	Left	44cm	52cm
NMA65A13/44R	Right	44cm	52cm






Cosmetic Soft Foam Cover
Modular - transfemoral prosthesis with conical 3-step-bore (Ø 70mm/Ø 60mm/Ø 32mm)

Part Number	Side	Calf Circ.	Knee Height
5 NMA65A14/44L	Left	44cm	52cm
NMA65A14/44R	Right	44cm	52cm

Cosmetic Soft Foam Cover
Transfemoral prostheses for Kinegen Stream knee N3A2000, N3A2100, N3A2200 - with conical pre-bore and with step-bore extra long (Ø 75mm/Ø 32mm)

Part Number	Side	Calf Circ.	Knee Height
6 NMA62A152/40L	Left	40cm	52cm
NMA62A152/40R	Right	40cm	52cm
NMA62A152/44L	Left	44cm	52cm
NMA62A152/44R	Right	44cm	52cm

Foam Cosmesis



Cosmetic Soft Foam Cover - Modular, Transtibial and Knee Disarticulation Prostheses				
Part Number	Length	Top	Bottom	Bore
1 NMA60A7	50cm	18.5cm	16cm	$\varnothing 30\text{mm}$

Cosmetic Soft Foam Cover - Modular and Transtibial Prostheses. Mouldable at 100-130°C				
Part Number	Length	Top	Bottom	Bore
2 NMA60A8	50cm	16cm	16cm	$\varnothing 30\text{mm}$
3 NMA60A9	50cm	16cm	16cm	$\varnothing 34\text{mm}$

Cosmetic Soft Foam Cover - Above Knee with Step Bore ($\varnothing 60\text{mm}/32\text{mm}$)			
Part Number	Side	Knee Height	Calf Circ.
4 NMA60A42/36L NMA60A42/36R	Left	42cm	36cm
	Right	42cm	36cm
NMA60A42/40L NMA60A42/40R	Left	42cm	40cm
	Right	42cm	40cm
NMA60A42/44L NMA60A42/44R	Left	42cm	44cm
	Right	42cm	44cm
NMA60A47/36L NMA60A47/36R	Left	47cm	36cm
	Right	47cm	36cm
NMA60A47/40L NMA60A47/40R	Left	47cm	40cm
	Right	47cm	40cm
NMA60A47/44L NMA60A47/44R	Left	47cm	44cm
	Right	47cm	44cm
NMA60A52/36L NMA60A52/36R	Left	52cm	36cm
	Right	52cm	36cm
NMA60A52/40L NMA60A52/40R	Left	52cm	40cm
	Right	52cm	40cm
NMA60A52/44L NMA60A52/44R	Left	52cm	44cm
	Right	52cm	44cm

Cosmetic Soft Foam Cover - Above Knee (same version as NMA60A47 but with extra long step-bore)			
Part Number	Side	Knee Height	Calf Circ.
5 NMA61A47/36L NMA61A47/36R	Left	47cm	36cm
	Right	47cm	36cm
NMA61A47/40L NMA61A47/40R	Left	47cm	40cm
	Right	47cm	40cm
NMA61A47/44L NMA61A47/44R	Left	47cm	44cm
	Right	47cm	44cm

Cosmetic Soft Foam Cover for Modular - Transfemoral Prosthesis (Uniprox)			
Part Number	Side	Calf Circ.	Size
6 NF042016	Left	40cm	Small
	Right	40cm	Small
NF042018	Left	44cm	Large
	Right	44cm	Large

7
Cosmesis

7
Cosmesis