

Structural and Socket Components - 30mm & 34mm Heavy Duty Titanium

30mm Trulife Components Weight Limit - 205kg (unless otherwise stated)

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.

30mm Heavy Duty Titanium Structural Components

Part Number	Description	Weight Limit
1 NAAAHTI210	4-Hole Male Pyramid Adapter	205kg
2 NAAAHTI217	Splined 4-Hole Male Rotatable Pyramid Adapter	205kg
3 NAAAHTI218	Splined 4-Hole Female Rotatable Receiver Adapter	205kg
4 NAAAHTI221	35mm Long Double Ended Adapter	204kg
5 NAAAHTI225	Clamp Adapter	204kg
6 NAAAHTI200	Bonded Pylon 420mm	205kg
7 NAAASMP	Seattle Mounting Plate - Includes Expulsion Valve Insert	205kg

This lamination plate is supplied with an expulsion valve as standard, but can be used without, if required.

Note: The weight limit is 168kg for the heavy duty aluminium components above for patients of K4 activity level.

34mm APC Components Weight Limit - 227kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.

34mm Heavy Duty Titanium Structural Components

Part Number	Description	Weight Limit
1 NHD-100C	Heavy Duty Titanium Tube Clamp	227kg
2 NHD-110LONG	Heavy Duty Titanium Receiver / Pylon Adapter - Long (82mm) Pylon	227kg
3 NHD-105C	Heavy Duty Titanium Tube Clamp with 5mm Offset	227kg
4 NHD-132D	Heavy Duty Titanium Dual Receiver Clamp	227kg
5 NHD-740L	Heavy Duty Titanium 4-Prong Lamination Pyramid	227kg

