

## Structural and Socket Components - Interchangeable Adapters



### Components Weight Limit - 100-166kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.

### Special Adapters

Part Number	Description	Weight	Build Height	Weight Limit
1	NTI-SL70-PF Sleeve System Rotatable Pyramid, Female	96g	21mm	166kg
2	NTI-SL70-PM Sleeve System Rotatable Pyramid, Male	54g	6mm	166kg
3	NTI-SL75-RF Sleeve System Rotatable Receiver, Female	91g	29mm	166kg
4	NTI-SL75-RM Sleeve System Rotatable Receiver, Male	62g	13mm	166kg
5	NTI-600-SLA-PM Sleeve System Pyramid, Male Dovetail	85g	7mm	100kg
6	NTI-610-SLA-TF Sleeve System Pyramid, Male Thread, Female Dovetail	77g	16mm	100kg
7	NTI-615-SLA-CF Sleeve System Tube Clamp, Female Dovetail	91g	48mm	100kg
8	NTI-620-SLA-FF Sleeve System Dual Female Dovetail	88g	32mm	100kg
9	NTI-630-SLA-RF Sleeve System Receiver, Female Dovetail	94g	34mm	100kg
10	NTI-640-SLA-PF Sleeve System Pyramid, Female Dovetail	94g	8mm	100kg
11	NTI-650-SLA-RM Sleeve System Receiver, Male Dovetail	99g	21mm	100kg
12	NTI-655-SLA-CM Sleeve System Tube Clamp, Male Dovetail	96g	35mm	100kg

**Note:** Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.