

QuickFit™ Lacers

The Click Medical QuickFit Lacers are suitable for both orthotic and prosthetic applications, with low-profile lace guides to quickly dial in fit. Using the easy grip Boa dial, users can easily achieve uniform compression with a one-handed twist to adjust fit. The kits are available in four ready-to-install sizes, which can then be easily replaced on devices without the need to remove the rivets or refabricate.



Features and Benefits

- Available in 4 kit sizes: X-Small, Small, Medium and Large
- Nylon-coated steel lace used for high strength and low friction
- Low-profile dial and lace guide design
- Double-cap press-fit steel rivets are easy to install using the QuickFit tool
- Lacers have an even-tension pull proximal to distal, so there is no need to manually redistribute tension

Ordering Information

Size	Order No.
X-Small*	NCMOK0180-130-05
Small*	NCMOK0200-130-05
Medium*	NCMOK0220-130-05
Large*	NCMOK0240-220-05
QuickFit Lacer Starter Kit	NCMOK0260-130-05

**Note: Lacers are sold individually. One pair of footwear requires two kits (one per shoe).*

Service Parts Ordering Information

Product	Order No.
QuickFit Setter Tool	NCMTL001
QuickFit Rivets 100 Pack	NCMSP1040
QuickFit XS/SM/MD Dial & Lace Replacement Kits	NCMRP0200-130-05
QuickFit LG Dial & Lace Replacement Kits	NCMRP0200-230-05
Bulk Lace Kit	NCML0K019

Parts Information

Size	Details
X-Small*	Comes with a low-profile dial, lace and 4 guides. Rivets sold separately. Designed for paediatric footwear size 35 EU/2 UK and smaller.
Small*	Comes with a low-profile dial, lace and 6 guides. Rivets sold separately. Designed for footwear size 41 EU/7 UK and smaller.
Medium*	Comes with a low-profile dial, lace and 8 guides. Rivets sold separately. Designed for footwear size 42 EU/8 UK and larger.
Large*	Designed for AFOs and above the ankle boots. Comes with a dial capable of more power, lace and 10 guides. Rivets sold separately.
QuickFit Lacer Starter Kit	Includes everything required to start using QuickFit Lacers: 10 QuickFit Lacers (2 x X-Small, 4 x Small, 4 x Medium), QuickFit Tool and 100 pack of rivets.
QuickFit Setter Tool	Recommended installation tool. Customised for use with QuickFit Rivets. Properly secures rivets and provides appropriate hole-punch for creating rivet holes when needed.
QuickFit Rivets 100 Pack	Available in packs of 100. Rivet pack is required to install any QuickFit Lacer Kit. Rivets are not included with QuickFit Lacer Kits.
QuickFit Dial & Lace Replacement Kits	Kits available in sizes X-Small, Small, Medium and Large - includes QuickFit Dial and QuickFit Lace.
Bulk Lace Kit	10m of QuickFit nylon-coated steel lace.
Information on bulk ordering QuickFit lacers and rivets available - please contact your local Product Manager.	

To order contact customerservices@steepergroup.com

1

2

3

4

5

6

7

8

9

10