

The Click Medical QuickFit Lacers are an easy to use fastening solution, allowing micro-adjustments throughout the day. They can be retrofit to footwear, or used for other light closure applications in prosthetics and orthotics to offer a range of benefits for those with upper limb absence, or limited mobility.

Features and Benefits

- **One-Hand, Micro-Adjustment**
Uniform compression is easily achieved with a one-handed twist of the easy grip dials, allowing patients with limited fine motor skills to adjust fit in seconds throughout the day
- **Pre-Laced System**
Kit choice and fabrication is made simple with four ready-to-install sizes
- **Easily Replaceable Components**
Easily replace the dial and lace, with no need to remove rivets or re-fabricate device
- **Proximal to Distal Tension Pull**
Even-tension pull proximal to distal means no need to manually redistribute tension as with traditional lacing. Nylon-coated steel lace is low friction and high strength
- **Rivet On**
With double-cap press-fit steel QuickFit rivets, these are retrofittable to other systems and easy to install with the QuickFit Tool, along with durable, low-profile lace guides



Ordering Information

Order No.	Size	Details
NCMOK0180-130-05	X-Small*	Supplied with a low-profile dial, lace and 4 guides. Rivets sold separately. Designed for paediatric footwear size 35 EU/2 UK and smaller.
NCMOK0200-130-05	Small*	Supplied with a low-profile dial, lace and 6 guides. Rivets sold separately. Designed for footwear size 41 EU/7 UK and smaller.
NCMOK0220-130-05	Medium*	Supplied with a low-profile dial, lace and 8 guides. Rivets sold separately. Designed for footwear size 42 EU/8 UK and larger.
NCMOK0240-230-05	Large*	Designed for AFOs and above the ankle boots. Comes with a dial capable of more power, lace and 10 guides. Rivets sold separately.
NCMOK0165-230	Quickfit Releasing Lacer	Adds opening functionality to the QuickFit Lacer by adding releasing guides to one side so that the system can fully open for easy donning. Rivets sold separately.
NCMSP1040	Quickfit Lacer Rivet Set (100 Pack)	Required to Install QuickFit Lacers to any device or shoe.

*Note: Lacers are sold individually. One pair of footwear requires two kits (one per shoe).

Service Parts Ordering Information

Order No.	Details
NCMRP0200-130-05	QuickFit® XS/S/M Dial and Lace
NCMRP0200-230-05	QuickFit® LG Dial and Lace
NCMLOK019	Low-Mid Power Steel Lace 10m

Please Note: Legacy spare parts are available to support older kits, please contact Customer Services or your Product Manager for more information.