

# QuickFit™ Strap and Buckle

The Click Medical QuickFit Strap and Buckle are an easy to use fastening solution, allowing for microadjustments throughout the day. Each 'click' of the BOA dial decreases strap length by 1mm increments, empowering patients to control comfort for the perfect fit everytime.

The strap is prefabricated and can be installed in minutes along with the QuickFit Buckle. The Buckle fastens through a self-locating magnet, providing additional strength and ease of use for those with limited manual dexterity.

## QuickFit Strap Features and Benefits

- Micro-Adjustable**  
Empower patients to quickly and easily control fit and comfort throughout the day. Each "click" of the dial decreases the strap length by 1mm for the perfect fit every time
- Improved Outcomes**  
Providing reliable adjustability to optimise fit improves the efficiency of a correctly fitting device
- Reliable Function**  
QuickFit Straps are prefabricated, durable, and can be installed in minutes

## QuickFit Buckle Features and Benefits

- High Durability**  
The QuickFit Buckle design prevents the effects of wear and tear experienced with hook and loop fastenings
- Self-Locating Buckle**  
Powerful magnetic buckle clicks itself into place and holds strong throughout the day, promoting ease of use for those with limited manual dexterity or visual impairment
- One-Handed Use**  
Opening and closing can be easily performed with one hand



# Structural and Socket Components - Click Medical

## Ordering Information

Product	Size (Width)	Order No.	Details
QuickFit Strap	3.8cm /1.5"	NCMOK1410-130-05	<ul style="list-style-type: none"><li>Hook and loop closure</li><li>10cm of adjustable travel</li><li>Usable strap length: Max = 75cm, Min = 20cm</li><li>Cut to size</li></ul>
QuickFit Buckle	3.8cm /1.5"	NCMOK1415-000-05	<ul style="list-style-type: none"><li>Easily added on to the QuickFit Strap</li><li>Adds 10cm to strap usable length</li><li>5mm profile</li><li>70kg holding capacity</li></ul>

