

1

QuickFit™ Straps & Buckles

2

The Click Medical QuickFit Straps and Buckles are an easy to use fastening solution, allowing for micro-adjustments throughout the day. Each 'click' of the dial decreases strap length by 1mm increments, empowering patients to control comfort for the perfect fit everytime.

3

The straps are prefabricated, durable, and can be installed in minutes along with the QuickFit buckle. The buckle fastens through a self-locating magnet, providing additional strength and ease of use for those with limited manual dexterity.

4



5

6

7

8

9

10



Features and Benefits

- Both components can be cut to size
- 70kg holding capacity
- QuickFit Buckle easily adds on to the QuickFit Straps, adding 10cm to the usable length
- Provides reliable adjustability to optimise fit, improving the efficiency of a correctly fitting device
- The QuickFit Buckle prevents the wear and tear experienced with hook and loop fastenings



Self-locating magnetic buckle

Ordering Information

Description	Size	Part Number
Quick Fit Strap	1" / 2.5cm	NCMOK1400-130-05
	1.5" / 3.8cm	NCMOK1410-130-05
	2" / 5cm	NCMOK1420-130-05
QuickFit Buckle (sold separately)	1" / 2.5cm	NCMOK1405-000-05
	1.5" / 3.8cm	NCMOK1415-000-05
	2" / 5cm	NCMOK1425-000-05