

RevoFit® Lamination Kit

The RevoFit® Lamination Kit allows you to create built-in user-controlled adjustability anywhere on a socket, letting users make incremental changes to their socket fit and volume throughout the day without the need to add or remove socks.

For use with all commonly used lamination materials, the RevoFit uses a tried and tested system of manufacture. This kit has been formulated to resist the repetitive stress found in an adjustable socket, incorporating a high durability (HD), fatigue-resistant lace. The system is well established and empowers users by allowing adjustability; it is complimented by the user-friendly Click Reel providing additional benefits including a balance of precise adjustment and power, offering up to 7:1 mechanical advantage.

Waterproof*

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Adjustable global and/or targeted compression
- Micro/macro adjustments to increase and to decrease socket volume
- Easy adjustments through clothing
- Optional grip rings to suit user dexterity needs
- Rapid donning & doffing for challenging residual limb shapes
- User Weight Limit: 135kg



Spend less time on remakes and refits.



Build adjustability, anywhere.

Ordering Information

RevoFit® Lamination Kit Part No.

NCMRF-200-07-01

Lamination Kit Includes:

- 1x Lamination Tube (1.8m)
- 1x Reel Tool
- 1x Lam. Collar



RevoFit® Additional Parts Information

Product	Part No.	Details
RevoFit® Lamination Kit with Click® Reel	NCMRF-200-07-10	10x Click® Reels, 10x Lamination Collars, 1x Spool of HD Lace (25m), 5x Metal Lace Feeders, 20x Plastic Lace Feeders, 3x Reel Tools, 10 Lamination Tubes (1.8m)
Click® Reel and HD Lace Kit	NCMRP-002-07	Click Reel, HD Lace (2m), Reel tool, 3x plastic lace feeders, 1x metal lace feeder. Note: This is not for replacing a BOA® dial. Please see NCMRP0100-320-05 for replacing a BOA® dial.
HD Lace Kit	NCMRP-002	HD Lace (2m), Reel tool, 3x plastic lace feeders, 1x metal lace feeder, red retainer tab
Lamination Collar	NCMRP-200	Replacement Lamination Collar for RevoFit & RevoLock Lamination Kits with the Click Reel only.
Lamination Tube (1.8m)	NCMRP-201	Replacement Lamination Tube (1.8m) for all RevoKits with the Click Reel and BOA® dials.
HD Lace Spool (25m)	NCMBL-002	HD (High Durability) Lace spool (25m). Compatible with all Kits including Click Reel and BOA dials.
Comfort Grip Ring (Black 5-Pack)	NCMAC-100-05	Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Especially suitable for patients with dexterity issues.
Power Grip Ring (Black 5-Pack)	NCMAC-101-05	Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Especially suitable for patients with dexterity and hand strength issues.

Please Note: Legacy spare parts are available to support older kits, please contact Customer Services or your Product Manager for more information.

Technical Information

RevoSurface® Tool Kit		Lace	2m HD (High Durability) Lace
Build Height of Click® Reel	14mm visible profile, 24mm total	User Weight Limit	300 lb/135kg
Reel Diameter	42mm		
Reel Weight	71g		
Collar Diameter	42mm		
Tube/Sheath Length	1.8m		
Tube Outer Diameter	3.2mm		



RevoFit Gap Socket

RevoFit Gap Socket

RevoFit 3D Printed Socket