

RevoLock® Align is a lanyard suspension system providing users with unrivaled control over the lanyard tension at all times, including when donning and doffing of a socket over a liner. The latest release introduces a number of additional features. These features include an optional Slider kit (to maintain proximal cable tension in sitting, even after distal tension is released), and a Quick Connect Snap insert (allows the lanyard to move more freely and the liner to be easily (dis) engaged). This kit is commonly used in both upper and lower limb devices, and allows the users alignment to be accommodated within the socket. RevoLock offers all the benefits associated with lanyard suspension.

The RevoLock® 4-Hole Kit offers an alternative way of applying the familiar RevoLock technology with industry standard 4 hole alignment systems.

The RevoLock Lanyard V2 Kit is designed to be similar to the original RevoLock Lanyard Kit, including the familiar BOA® dial and the former distal housing/threaded insert. The RevoLock Lanyard V2 restores a well-established adjustable suspension solution that meets with user expectations in a simple, easy to follow way.

Waterproof*

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

Align Kit

- Compact, suitable for users with residual limbs of a variety of shapes and of the usual sizes
- Transferable alignment
- Locked and hidden chord
- 38mm build height
- Reduced pistoning
- Eliminates locking pin alignment frustrations
- Slider kit allows user to lock the Click Reel power mode on demand
- No tool form required for manufacture
- User Weight Limit: 135kg

4-Hole Kit

- Industry standard 4-hole interface
- Locked and hidden chord
- 25mm build height
- Reduced pistoning
- Eliminates locking pin alignment frustrations
- Slider kit allows user to lock the Click Reel power mode on demand
- RevoLock® Tool form required for fabrication, sold separately
- User Weight Limit: 135kg

Lanyard V2 with BOA Dial

- For transferable alignment
- Single or double-stage lamination
- Low 25mm build height
- Reduced pistoning
- Upper or Lower extremity sockets
- Lightweight system
- No tool form required

Ordering Information

RevoLock® Align Kit Part No.

NCMRL-301-07-01

Align Kit Includes:

- 2x RevoLock® LE Quick Connect Inerts & 1x Snap
- 1x RevoLock® Slider Lamination Kit
- 1x Align Housing & Fabrication Post
- 1x Lamination Tube (0.5m)
- 1x Reel Tool,
- 1x Lamination Collar
- 1x Click Reel,
- 1x HC Lace (2m)
- 3x Plastic Lace Feeders
- 1x Metal Lace Feeder

RevoLock® 4-Hole Kit Part No.

NCMRL-300-07-01

4-Hole Kit Includes:

- 2x RevoLock® LE Quick Connect Inerts & 1x Snap
- 1x RevoLock® Slider Lamination Kit
- 1x 4-Hole Puck
- 1x Lamination Tube (0.5m)
- 1x Reel Tool
- 1x Lam. Collar
- 1x Click® Reel
- 1x HC Lace (2m)
- 1x Tool Form Plug
- 3x Plastic Lace Feeders
- 1x Metal Lace Feeder

RevoLock Lanyard V2 with BOA Dial

NCMRL3000-220

RevoLock® Lanyard V2 with BOA® Dial Replacement Kits

NCMRP-310	RevoLock V2 BOA Dial and HC Lace Kit (Port C)
NCMRP-310A	RevoLock V2 BOA Dial and HC Lace Kit (Port A)
NCMRP-312	RevoLock V2 HC Lace Kit

RevoLock® Tool Form Part No.

NCMRL-310-00-01

Please Note: Tool form is required and must be ordered in addition to the kit

Technical Information

	RevoLock® Align Kit	RevoLock® 4-Hole Kit	RevoLock® 4-Hole Kit
Build Height of Click® Reel	14mm visible profile, 24mm total	14mm visible profile, 24mm total	19mm
Reel Diameter	42mm	42mm	36mm
Reel Weight	71g	71g	14g
Housing Diameter	45mm (Distal)	45mm	30mm
Housing Height	35mm (Distal)	24.5mm	25mm
Housing Weight	28g (Distal)	116g	21.8g
Snap/Threaded Insert Weight	33g	33g	5.9g
Tube/Sheath Length	0.5m	0.5m	0.5m
Lace	2m HC (High Capacity) Lace	2m HC (High Capacity) Lace	2m HC (High Capacity) Lace
User Weight Limit	300 lb/135kg	300 lb/135kg	200 lb/100kg



Components

Components