


Structural and Socket Components - Special Adapters

Components Weight Limit - 100-166kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.



Special Adapters

Part Number	Description	Weight Limit
1 NMA11A15/A	Aluminium Tube Clamp Adapter with Pyramid Ø 30mm	100kg
2 NMA10A41/T	Titanium Female Lamination Adapter with Shift	100kg
3 NMA10A43/T60	60mm Titanium Double Adapter with Shift	100kg
NMA10A43/T75	75mm Titanium Double Adapter with Shift	100kg
4 NMA10A45/T	Titanium Double Pyramid with Shift	100kg
5 NMA11A14/A	Aluminium Double Adapter with Pyramid	100kg
6 NTI-125P	Titanium Pyramid and Receiver	166kg
7 NMA10A60/A41	Aluminium Double Pyramid Adapter 15mm	100kg
NMA10A60/A50	Aluminium Double Pyramid Adapter 23mm	100kg
8 NTI-115DM	Titanium Dual Pyramid Adapter 16mm	166kg
9 NMA11A16/T	Titanium Pyramid Adapter with 30mm Core	125kg
NMA11A16/S	Stainless Steel Pyramid Adapter with 30mm Core	150kg
10 NTI-112P	Titanium 30mm Pylon with Pyramid, Short	166kg
11 NMA10A89/T	Titanium Four Hole Female Adapter Offset	100kg

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.