

## Structural and Socket Components - Waterproof with Drainage

### Weight Limit - 166kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight limit of the build.



### Aluminium 30mm Structural Components

Part Number	Description	Weight Limit
1	NTI-100C-DH Titanium Tube Clamp With drain holes and contains water resistant fasteners.	166kg
2	NTI-10TC2-DH 262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
3	NTI-17TC2-DH 440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
4	NTI-10TC2-DH-BLACK 262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
5	NTI-17TC2-DH-BLACK 440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
6	NTI-740L 4 Prong Titanium Laminating Adapter with Pyramid	166kg

**Note:** All components listed above are not only waterproof, but have drainage slots to ensure any water ingress into the tube can escape. This makes these components ideal for use on shower limbs, or for any application wherever there is submersion in, or frequent contact with water.