

Structural and Socket Components - Waterproof with Drainage

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.



Titanium 30mm Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	NTI-10TC2-DH-SILVER	APC	262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
2	NTI-17TC2-DH-SILVER	APC	440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
3	NTI-10TC2-DH-BLACK	APC	262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
4	NTI-17TC2-DH-BLACK	APC	440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
5	NTI-100C-DH	APC	Titanium Tube Clamp With drain holes and contains water resistant fasteners.	
6	NTI-740L	APC	4 Prong Titanium Laminating Adapter with Pyramid	166kg
7	NTI-740-S	APC	4 Prong with Pryamid, Slots	166Kg

Note: The drainage slots included ensure these components are ideal for use on shower limbs, or for any application wherever there is submersion in, or frequent contact with water.

Note: See pages 246-247 for the Mix and Match APC Alignment Adapters Component Build Height Information

Manufacturer: APC = [American Prosthetic Components](#)