

# Breeze

Breeze provides K2 users with a cost effective, yet highly functional foot with a variety of features utilising College Park's high-end technology. The water friendly system has an optional water drainage port, corrosion-resistant materials and special modifications to clear water from the entire prosthetic system.

**CES CP F 11:** Use this code to read the Breeze clinical evaluation summary on the Steeper website.

## Features and Benefits

- Low profile
- IntelliWeave composite for assisted movement
- Toe roll-over assist Enviroshell
- Removable sandal toe foot shell with antimicrobial properties
- Water friendly with optional drainage port
- Integrated pyramid



Waterproof\*



## Ordering Information

Trial Period - This foot has a 30 day trial period.

### Impact Level Descriptions

**Low:** Daily activities include mostly level ground walking, moving around in the home and the community.

Endo IP - NBE ■●EN▲ Integrated Pyramid

- - L or R for Left or Right
- - Colour of Cosmesis (C = Caucasian, T = Tan, B = Brown)
- ▲ - Foot Size in cm (21-30)

Includes foot, removable foot shell and CPI sock. See order forms on page 36 for additional ordering information.

## Product Specification

Mounting	Weight Limit	Sizes	Impact Level	Build Height (cm)	Weight	Warranty	Foot Shell
Endo (IP)	100kg (21-25cm) 125kg (26-30cm)	21-30cm	L M H	6.2 (21-25cm) 6.5 (26-30cm)	523g (Size 25cm including foot shell)	2 years Foot shell - 3 months	C T B Caucasian, Tan or Brown Standard width

\*Waterproof note: Rinse with fresh water and dry following exposure to salt water or chlorinated water.

## Shelltread - Ordering Information

Shelltread	Size
NBET ■●▲	21-30cm

Spare Parts Ordering Guide

- = "L" or "R" for Left or Right
- = Colour (N = Neutral, J = Jet Black)
- ▲ = Foot size in cm (21-30)

# Feet - College Park - Breeze

## Technical Information

