

The Orbit foot was designed to provide smooth rotation and rollover. College Park's unique, built-in torsion feature offers a total rotation range of 40° within an incredibly low profile. The Orbit incorporates College Park's composite expertise in a carbon fibre design to provide efficient energy return, which can improve comfort during sports, including golf or tennis.

Waterproof\*

\*Waterproof note: Following exposure to chlorinated and salt water, rinse immediately with fresh water and dry.

Features and Benefits

- Adjustable torsion with integrated adapter giving a total rotation range of +/- 20°
- Effective energy return
- Low profile design
- Strong carbon fibre base provides a dynamic response during low to high impact activities
- Waterproof\* (fresh water only)



Smooth roll over



+/- 20° Adjustable torsion



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

**Moderate:** Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

**High:** Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Ordering Example - NOR-LCEN-24-G3						Weight (kg)	Low/Mod	High
Part ID	Side	Shell Colour	Mounting	Size	Firmness Category			
NOR	L		EN		G	0-65	1	2
NOR	R		EN		G	65-82	2	3
						82-100	3	4
						101-125	4	5
						126-150	5	6

  

Caucasian	C
Tan	T
Brown	B
Jet Black	J

21-30cm

Spare Parts Ordering Information

Foot Shell	Size (cm)	Wedge Kit	Size (cm)
NHZ-S ■ ● ▲	21-30	NHZ-HWK-2125	21-25
		NHZ-HWK-2630	26-30

  

CAP	Size (cm)	CPI Spectra Sock	Size (cm)
NBE-C ■ ● -2125	21-25	NCPI-SS-SM	Small (19-25cm)
NBE-C ■ ● -2630	26-30	NCPI-SS-LG	Large (26-31cm)
		NCPI-SSK-SM	Small 6 pack (19-25cm)
		NCPI-SSK-LG	Large 6 pack (26-31cm)

  

Shelltread	Size (cm)
NHZ-T ■ N-▲	21-30
Neutral	

■ = "L" or "R" for Left or Right  
 ● = Colour (C = Caucasian, T = Tan, B = Brown, J = Jet Black)  
 ▲ = Foot size in cm (21-30)

Product Specification

Mounting	Sizes	Build Height	Weight Limit
Endo (IP)	21-30cm	21-25cm 9.4 cm 26-30cm 10 cm	150kg

  

Weight of Unit	Impact Level	Warranty	Foot Shell
895g Size 26 inc. foot shell	M H	3 years foot shell - 6 months	C T B J Caucasian, Tan, Brown or Jet Black Standard Width

