

Delivering a lower profile than the Soleus, the upgraded Velocity features coupled toe springs that work together to provide a similarly progressive and smooth roll-over.

Velocity is a high-performance, low maintenance foot that is simple to align and set up. The additional inserts provided allow fine tuning of the foot function.

CES CPI F10: Use this code to read the Velocity clinical evaluation summary on the Steeper website.

Features and Benefits

- Low-profile Integrated Spring Technology (IST)
- Multi-axial function
- Intelliweave® composite to bring maximum strength and flexibility
- Easily adjustable tuning with three different heel wedges and a toe wedge provided



Wedges for customisation



Integrated spring technology

Cosmetics



CAP



Solid toe foot shell



Ordering Information

These devices do not have a trial period as standard. If you wish to discuss the Velocity please contact your local product manager for advice prior to placing your order.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Ordering Example: NVL-LCMEN-29N-G2

Part ID	Side	Shell Colour	Housing/ Mounting	Size	Width	Firmness Category
NVL	L		MEN			G
NVL	R		MEN			G

Includes CPI Spectra Sock and foot shell.


Caucasian	C	Brown	B	25-30cm	Narrow	N	1-6
Tan	T	Jet Black	J		Wide	W	See Chart Below

Firmness Category

Weight (kg)	Sizes 25-26cm		Sizes 27-30cm	
	Low - Mod Impact	High Impact	Low - Mod Impact	High Impact
0-45	1	1	1	1
46-54	1	2	1	2
55-63	1	2	1	2
64-72	2	3	1	2
73-81	2	3	2	3
82-91	3	4	2	3
92-100	4	5	3	4
101-113	5	6	4	5
114-125	-	-	5	6

Spare Parts Ordering Information

Foot Shell	Size (cm)
NV-S	25-30



CAP	Size (cm)
NV-C	25-26
NV-C	27-30



Shelltread	Size (cm)
NV-T	25-30




- = "L" or "R" for Left or Right
- = Colour (C = Caucasian, T = Tan, B = Brown, J = Jet Black)
- ▲ = Foot size in cm (25-30)
- ◆ = N = Narrow, W = Wide

Wedge Kit	Size (cm)
NVL-HWK-2125	21-25
NVL-HWK-2630	26-30




CPI Spectra Sock	Size (cm)
NCPI-SS-SM	Small (19-25cm)
NCPI-SS-LG	Large (26-31cm)
NCPI-SSK-SM	Small 6 pack (19-25cm)
NCPI-SSK-LG	Large 6 pack (26-31cm)



Product Specification

Mounting	Sizes	Weight Limit	Build Height
Endo (IP)	25-26cm	113kg	12cm
	27-30cm	125kg	13.2cm
Weight of Unit	Impact Level	Warranty	Foot Shell
631g Size 26 inc. foot shell		3 years Foot shell - 6 months	C T B J Caucasian, Tan, Brown or Jet Black



To order contact customerservices@steepergroup.com