

This light-weight and compliant foot has been designed for users with low to moderate levels of impact.

This foot has a low profile, similar to the Trés and Tribute feet. It also has built-in springs for a smooth roll over and added security for users, especially on uneven surfaces. Its action is softer than the Tribute and Trés, making it ideal for those patients who need a foot that protects them from excessive socket interface forces whilst also improving stability.

This durable yet lightweight foot is packaged within a low maintenance design, complete with a removable foot shell.

CES CPI F09: Use this code to read the Celsus clinical evaluation summary on the Steeper website.

Features and Benefits

- Integrated spring system and Intelliweave® technology provides excellent multi-axial movement with a natural centre of rotation; providing ground compliance and increased stability
- Low profile
- Removable foot shell with sandal toe
- Custom gait matching with pre-installed heel wedge and fulcrum pad to allow precise levels of resistance for maximum comfort



Multi-axial response

Cosmetics



Removable foot shell with sandal toe



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around the home and the community.

Ordering Example: NCS-LCEN-21-G2

Part ID	Side	Shell Colour	Mounting	Size	Firmness Category
NCS	L	<input type="text"/>	EN	<input type="text"/>	G <input type="text"/>
NCS	R	<input type="text"/>	EN	<input type="text"/>	G <input type="text"/>

Includes foot shell and CPI Spectra Sock.

Caucasian	C
Tan	T
Brown	B

21-30cm

1-4

Firmness Category				
Weight (kg)	0-63	64-81	82-100	101-136
Size (cm)	21-30		25-30	
Low Impact	1	2	3	4

Spare Parts Ordering Information

Foot Shell	Size (cm)
NCS-S ■ ● ▲	21-30

Shelltread	Size (cm)
NCS-T ■ N-▲	21-30cm

■ = "L" or "R" for Left or Right
 ● = Colour (C = Caucasian, T = Tan, B = Brown)
 ▲ = Foot size in cm (21-30)

Neutral

CPI Spectra Sock	Size (cm)
NV-SS-2124	21-24
NV-SS-2526	25-26
NV-SS-2730	27-30
NV-SSK-2124	21-24 (6 pack)
NV-SSK-2526	25-26 (6 pack)
NV-SSK-2730	27-30 (6 pack)

Product Specification

Mounting	Size	Build Height	Weight Limit
Endo	21-24cm	5.7cm	100kg
	25-26cm	6.4cm	136kg
	27-30cm	6.6cm	

Weight	Impact Level	Warranty	Foot Shell
526g Size 26cm inc. foot shell		3 years Foot shell - 6 months	C T B Caucasian, Tan or Brown Standard width

Technical Information

