

Defender Paediatric

The Defender is a specifically designed foot to endure the physical demands of the young prosthesis user. The carbon fibre lay-ups are optimised to provide flexibility in a design that is robust and virtually maintenance-free so that users can feel confident over a range of high impact activities.

Features and Benefits

- Patented technology and eliminated bolt holes result in optimal strength and performance
- Uniquely optimised for the varied gait of a child
- Virtually maintenance-free design
- Available in 6 stiffness categories
- Provides the user with the ability to seamlessly transition from the paediatric Defender foot to an adult Highlander foot



Water resistant - safe for contact with water and submersion*

Cosmetics



Mid profile foot shell



Ordering Information

There is no trial period on this foot.

Ordering Example - NFS34P19L

Part Number	Category	Size	Left / Right
NFS3	4P	19	L


Includes foot shell and Spectra® Sock

Category Selection

Weight (kg)	Category
9-13	1P
14-18	2P
19-22	3P
23-27	4P
28-31	5P
32-36	6P
37-45	1A

Note: Category 1A allows a transition from the Defender Paediatric to the adult Highlander foot from size 22.

Technical Information



Product Specification

Build Height and Weight

Size (cm)	Build Height (cm)	Foot inc. Shell
19-21	7.6	225g

Weight Limit	Impact Level	Sizes (cm)	Warranty
45kg	MH	19-21	18 months Foot shell - 6 months

*Water note: Rinse with fresh water and dry following exposure to water, particularly salt or chlorinated.

Spare Parts Ordering Information



Foot Shell	Size
NFTC-2M-1 ▲ ● 4-R ■	19-21

Foot Shell Ordering Guide

- ▲ = Size in cm (19-21)
- = Skin Tone (L = Light, M = Medium, D = Dark)
- = "L" or "R" for Left or Right

Bumper Kit	Size (cm)
NKIT-00-12201-00	19-21

Spectra® Sock	Size (cm)	Sock Size (cm)
NSO-NPS-X0025-00	19-21	25

X denotes colour. Insert 1 for white, or 2 for black. **Note:** Default colour is black