

	Name	Impact Level	K Level	Features and Benefits	Foot Size (cm)	Build Height on size 25cm	Weight based on size 25cm inc foot shell	Weight Limit	Warranty	Part Number	Trial Period	Sandal Toe Available
Everyday	Senator		K2	<ul style="list-style-type: none"> <li>Lightweight, energy-return design</li> <li>Cost effective</li> <li>Low maintenance, durable design</li> <li>Available in 6 stiffness categories</li> </ul>	22-31	13cm	520g	136kg	2 Years	NVS1 ●▲■	None	No
	Sierra		K3	<ul style="list-style-type: none"> <li>Full length heel ensures exceptionally smooth gait</li> <li>Split keel design</li> <li>Integral male pyramid</li> <li>Available in 9 stiffness categories</li> </ul>	22-31	14.6cm	515g	166kg	3 Years	NFS1 ●▲■	30 day trial (only sizes 23-28cm, Categories 2-7)	Yes*
	Highlander		K2,3	<ul style="list-style-type: none"> <li>Low build height with integral male pyramid</li> <li>Split keel design for excellent inversion/eversion</li> <li>Increased weight limit</li> </ul>	22-31	10.8cm	470g	166kg	3 Years	NFS3 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	Highlander Max		K3		26-31	12.1cm (26cm)	590g (26cm)	226kg	1 Year	NFS3 (Categories 10-12)	Highlander Max - None	
	DynAdapt		K3,4	<ul style="list-style-type: none"> <li>Slim profile, carbon fibre foot system with a slender, anatomical design</li> <li>Multi-axial function for maximum comfort</li> <li>Full length split heel spring provides smooth and compliant transition from heel strike to foot flat</li> <li>Full length split keel gives an uninterrupted rollover with good compliance and a dynamic response through to toe off</li> </ul>	22-31	15.8cm	535g	166kg	3 Years	NF10 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	Renegade AT		K3,4	<ul style="list-style-type: none"> <li>Carbon fibre Z-shaped heel is durable and maintenance free</li> <li>Patented Z-Shock Technology™</li> <li>Lightweight solution at 1/3 the weight of traditional shock foot products</li> <li>EnduraCore Technology to improve durability</li> <li>Split-toe design delivers multi-axial function with +/- 15 degrees of inversion/eversion</li> </ul>	22-31	16.2cm	625g	166kg	3 Years	NR11 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	Renegade AT LP				12.7cm	625g	NR16 ●▲■					
	Maverick Comfort AT		K2,3	<ul style="list-style-type: none"> <li>Low build height for those with clearance challenges</li> <li>Split keel for increased ground compliance and versatility</li> <li>Increased durability without compromising the desired energy storage/return</li> <li>Smooth roll-over during stance phase of gait</li> <li>Lightweight fibreglass</li> </ul>	23-30	11.6cm	626g	166kg	3 Years	NF22 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
Maverick Xtreme		K3,4	<ul style="list-style-type: none"> <li>High impact, extreme user</li> <li>Regular or Sandal toe option</li> <li>High durability, flexibility and natural feel</li> <li>Lightweight fibreglass option</li> </ul>	23-30	16.2cm	650g	166kg	3 Years	NF11 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*	
Maverick Xtreme AT		K3,4	<ul style="list-style-type: none"> <li>High impact, extreme user</li> <li>Regular or Sandal toe option</li> <li>High durability, flexibility and natural feel</li> <li>Lightweight fibreglass option</li> <li>Full length split keel for excellent inversion/eversion on uneven terrains</li> </ul>	23-30	16.2cm	646g	166kg	3 Years	NF21 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*	
Speciality	Pacifica LP		K2,3	<ul style="list-style-type: none"> <li>Ultra low build height</li> <li>Integral male pyramid</li> <li>Split keel design</li> <li>Foot lies within the foot shell to improve cosmetic finishing</li> </ul>	22-31	4.7cm	465g	166kg	3 Years	NFS4 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	LP Symes		K2,3	<ul style="list-style-type: none"> <li>Carbon fibre for flexibility and comfort</li> <li>Designed specifically for the Symes amputee</li> <li>Lightweight</li> <li>Direct lamination cup</li> </ul>	22-31	4.3cm	420g	147kg	3 Years	NLP2 ●▲■	None	No
	Runway		K2,3	<ul style="list-style-type: none"> <li>Push button or Hex Key adjustment</li> <li>Anatomical 'gliding' ankle design ensures that alignment remains consistent, for smooth gait at all heel heights</li> <li>Carbon fibre design provides energy return</li> <li>Safe for occasional contact with water, but not intentional frequent contact or submersion</li> </ul>	22-28	9.8cm	690g	116kg	18 Months	NRS2 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Standard
	Thrive		K3,4	<ul style="list-style-type: none"> <li>Dual keel design to ensure consistent response under loading</li> <li>Full length primary keel selected against the patient's weight and activity level</li> <li>Split primary keel provides good inversion/eversion ground compliance</li> <li>Good shock absorption at heel strike and mid-stance</li> </ul>	22-31	14.5cm	630g	166kg	3 Years	NFS5 ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	Freestyle Swim Foot		K2,3,4	<ul style="list-style-type: none"> <li>Users can now walk comfortably to and from the swim area</li> <li>Adjust to 70° plantar flexion for swimming with the simple push of a button</li> <li>Utilise a fin to gain swimming speed and decrease energy expenditure</li> <li>The extended foot position balances the swimming action and makes it more efficient</li> </ul>	22-31	10.2cm	680g	116kg	1 Year	NLP2W2 ●▲■	None	No
Extreme Sports	Nitro Running Foot		K3,4	<ul style="list-style-type: none"> <li>Individually hand crafted ensuring optimal vertical displacement and energy return</li> <li>Integral male pyramid for easy alignment and manufacture</li> <li>Compatible with exoskeletal design</li> <li>Available in short or tall profile</li> <li>Stiffness tailored to amputee weight and running experience</li> </ul>	N/A	24.1cm (tall)	300g	166kg	1 Year	NFS6 ●	None	N/A
	Catapult		K3,4	<ul style="list-style-type: none"> <li>Dual spring technology increases the amount of energy stored and optimises the timing of its release</li> <li>Energy propulsion allows users to expend less energy in a more stable and balanced design</li> <li>Pyramid attachment facilitates alignment changes and the use of a single prosthetic socket for all activities</li> </ul>	N/A	24.1cm	510g	166kg	1 Year	NFX6 ●	Trial available on request	N/A
	Slalom Ski Foot		K3,4	<ul style="list-style-type: none"> <li>Eliminates the need for a ski boot</li> <li>Provides necessary flex and shock absorption for fluid motion</li> <li>Better control and range of motion increases endurance, allowing more time on the slopes</li> <li>Slip resistant rubber sole provides stability in snow and icy conditions</li> <li>Water resistant</li> </ul>	Custom Made	14cm	1055g	114kg	1 Year	To order contact <a href="mailto:customerservices@steepergroup.com">customerservices@steepergroup.com</a>	None	N/A
Hydraulic Ankles	Kinterra		K2,3	<ul style="list-style-type: none"> <li>Articulating ankle delivers reduced socket pressures for increased comfort</li> <li>Increased ground contact and a more symmetrical gait</li> <li>Improved stability and safety, especially on slopes and varied terrain</li> <li>Dorsi-Assist spring aids ground clearance in swing phase</li> <li>Safe for occasional contact with water, but not intentional frequent contact or submersion</li> </ul>	22-31	12cm	795g	125kg	3 Years	NROM ●▲■	30 day trial (only sizes 23-28cm, categories 2-7)	Yes*
	Kinnex		K2,3	<ul style="list-style-type: none"> <li>Microprocessor controlled resistance</li> <li>Dynamic lock</li> <li>Rapid Response™ valve technology</li> <li>Apple/Android compatible user app</li> <li>30 degrees of ankle articulation</li> <li>Waterproof and IP67 rated</li> </ul>	24-30	18.4cm	1432g	125kg	3 - 5 Years	NF14 ●▲●	● - insert 'RR' for regular toe right, 'RL' for regular toe left, 'SR' for sandal toe right, 'SL' for sandal toe left	Trial available on request
Paediatric	Defender		K2,3,4	<ul style="list-style-type: none"> <li>Patented technology and eliminated bolt holes result in optimal strength and performance</li> <li>Uniquely optimised for the varied gait of a child</li> <li>Virtually maintenance-free design</li> <li>Available in 6 stiffness categories</li> <li>Provides the user with the ability to seamlessly transition from the paediatric Defender foot to an adult Highlander foot</li> </ul>	16-21	6.1cm (Size 16-18) 7.6cm (Size 19-21)	175g (Size 16-18) 225g (Size 19-21)	45kg	18 Months	NFS3 ●▲●	None	No

= Water resistant (Designed to be water resistant, safe for contact with water and submersion)

= Waterproof (Designed for swimming and frequent use in water)

**Order Guide**

● = Category (1-9)\*

▲ = Size in cm (22-31)

■ = "L" or "R" for Left or Right

\*Sandal toe option only available in sizes 22-28.

Note: For sandal toe option add an "S" before the side of foot e.g when ordering Kinterra category 1, size 22 sandal toe left foot, the part number would be **NROM122SL**

Note: Category range differs for each foot, please refer to the Steeper catalogue for category tables.

Note: After use in water, the feet should be rinsed/cleaned.