

Freedom Highlander® MAX



Feet - PROTEOR - Freedom Highlander MAX

The Highlander Max enables the benefits of the Highlander to be made available to amputees over 166kg. With an industry highest weight limit of 226kg, this foot can be used in conjunction with the Heavy Duty Component Kit (NAC 202) to provide a prosthesis to virtually any transtibial amputee.

Water-resistant*

***Water-resistant:** Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.


Features and Benefits

- Low build height with integral male pyramid
- Split keel design for excellent inversion/eversion
- Available up to weight limit 226kg
- Available in 3 stiffness categories

Cosmetics



Mid profile foot shell



Sandal toe: sizes 26-28cm only



Ordering Information

Note: Highlander Max is a special order item, therefore please allow a 10-14 day lead time. There is no trial on this product.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

Ordering Example - NFS31126L

Part Number	Category	Size (26-31)	Left/Right
NFS3	11	26	L

Sandal Toe Ordering Example - NFS31126SL

Part Number	Category	Size	Sandal Toe	Left/Right
NFS3	11	26	S	L

Includes foot shell and Spectra® Sock

Category Selection - Highlander MAX (kg)

166-190kg	10
191-217kg	11
218-226kg	12



34mm Heavy Duty Component Kit (Max Only)	Weight Limit
NAC-202	226kg

Two 34mm Tube Clamps
4-Prong Connector, 34mm Pylon

Replacement - 4-Prong Connector

NACC-00-20202-00

Product Specification

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
226kg	L M H	590g Size 26cm	26-31	12.1 Size 26cm	1 year Foot shell - 6 months

Spare Parts Ordering Information

Part Number	Foot Shell Type
NFTC-3M-0 ▲●4-R ■	Standard, with CAP
NFTC-2M-1 ▲●4-S ■	Sandal Toe, without CAP



Cosmetic Attachment Plate

NCAP-3M-PO- ▲●4- ■



Spare Parts Ordering Guide

- ▲ = Size in cm (22-31 = Regular, 22-28 = Sandal Toe)
- = Skin Tone (L = Light, M = Medium, D = Dark)
- = "L" or "R" for Left or Right


Spectra® Sock	Colour	Size (cm)	Sock Size (cm)
NSO-NPS-20030-00	B	22-25	30
NSO-NPS-20035-00	B	26-28	35
NSO-NPS-20040-00	B	29-31	40

Colours: B = Black

Bumper Kit	Size (cm)
NKIT-00-11304-00	26-28
NKIT-00-11305-00	29-31

Build Height and Weight

Size (cm)	Build Height (cm)	Foot inc. Shell
26-28	12.1	590g
29-31	30	750g



Technical Information

Build Height and Weight		
Size (cm)	Build Height (cm)	Foot inc. Shell
22-25	10.8	470g
26-28	12.1	590g
29-31	30	750g