

# Freedom Pacifica LP®



The Pacifica LP offers proven performance in an ultra-low profile design, which requires less than 50mm of clearance for most commonly used sizes. The split-keel design provides stability-enhancing inversion and eversion as well as being lightweight and durable.



**\*Water-resistant:** Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.


## Features and Benefits

- Ultra low build height
- Integral male pyramid
- Split keel design
- Foot lies within the foot shell to improve cosmetic finishing
- Available in 9 stiffness categories

## Cosmetics



Foot shell



Sandal toe: sizes 22-28 only



# Feet - PROTEOR - Freedom Pacifica LP

## Ordering Information

**Trial Period** - This foot has a 30 day trial period on sizes 23-30 and categories 2-8 only.

## Impact Level Descriptions

**Low:** Daily activities include mostly level ground walking, moving around in the home and the community.

**Moderate:** Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

**High:** Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

**Note:** Although the 'high' impact level states 'running' as an appropriate activity, the Pacifica LP is not ideal for running.

Ordering Example - NFS4526L			
Part Number	Category	Size (22-31)	Left/Right
NFS4	5	26	L

Sandal Toe Ordering Example - NFS4526SL				
Part Number	Category	Size	Sandal Toe	Left/Right
NFS4	5	26	S	L

Includes foot shell and Spectra® Sock

Category Selection			
Weight (kg)	Impact Level		
	Low	Mod	High
44-52	1	1	2
53-59	1	2	3
60-68	2	3	4
69-77	3	4	5
78-88	4	5	6
89-100	5	6	7
101-116	6	7	8
117-130	7	8	9
131-147	8	9	-
148-166	9	-	-

## Product Specification

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
166kg		465g Size 25cm	22-31	4.7 Size 25cm	3 years Foot shell - 6 months

## Spare Parts Ordering Information


Part Number	Foot Shell Type
NFTC-3M-0 ▲●4-R ■	Standard, with CAP
NFTC-2M-1 ▲●4-S ■	Sandal Toe, without CAP
<b>Cosmetic Attachment Plate</b>	
NCAP-3M-PO- ▲● 4-■	

**Spare Parts Ordering Guide**  
 ▲ = Size in cm (22-31 = Regular, 22-28 = Sandal Toe)  
 ● = Skin Tone (L = Light, M = Medium, D = Dark)  
 ■ = "L" or "R" for Left or Right

Spectra® Sock	Colour	Size (cm)	Sock Size (cm)
NSO-NPS-20030-00	B	22-25	30
NSO-NPS-20035-00	B	26-28	35
NSO-NPS-20040-00	B	29-31	40

Colours: B = Black

Bumper Kit	Size (cm)
NKIT-00-11323-00	22-25
NKIT-00-11324-00	26-28
NKIT-00-11325-00	29-31



## Technical Information

Build Height and Weight		
Size (cm)	Build Height (cm)	Foot inc. Shell
22-25	4.7	465g
26-28	5	570g
29-31	6	760g